

Product Datasheet

HRAS Antibody - BSA Free

NBP2-42864

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-42864

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-42864



NBP2-42864

HRAS Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS (pH 7), 20% Glycerol
Target Molecular Weight	21 kDa

Product Description


Description	Novus Biologicals Rabbit HRAS Antibody - BSA Free (NBP2-42864) is a polyclonal antibody validated for use in IHC, WB, ICC/IF and IP. Anti-HRAS Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3265
Gene Symbol	HRAS
Species	Human, Mouse, Rat
Reactivity Notes	Xenopus laevis (100%).
Specificity/Sensitivity	This antibody is specific for H-Ras protein, and it does not cross react with human N-Ras or K-Ras protein.
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human HRAS. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:10000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunoprecipitation Assay dependent, Immunohistochemistry-Paraffin 1:100-1:1000

Images


Western Blot: HRAS Antibody [NBP2-42864] - Various whole cell extracts (30 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with H-RAS antibody diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.




Immunocytochemistry/Immunofluorescence: HRAS Antibody [NBP2-42864] - Analysis of paraformaldehyde-fixed HeLa, using H-Ras antibody at 1:200 dilution.




Immunohistochemistry-Paraffin: HRAS Antibody [NBP2-42864] - Mouse intestine. H-RAS antibody dilution: 1:500. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min.




Western Blot: HRAS Antibody [NBP2-42864] - Non-transfected (-) and transfected (+) 293T whole cell extracts (30 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with H-Ras antibody diluted at 1:5000.




Immunohistochemistry-Paraffin: HRAS Antibody [NBP2-42864] - Mouse intestine. H-RAS antibody dilution: 1:500. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min.




Western Blot: HRAS Antibody [NBP2-42864] -Whole cell extract (30 ug) was separated by 12% SDS-PAGE, and the membrane was blotted with H-Ras antibody diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.



Western Blot: HRAS Antibody [NBP2-42864] -Various tissue extracts (50 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with H-Ras antibody diluted at 1:2000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.



Western Blot: HRAS Antibody [NBP2-42864] -Various whole cell extracts (30 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with H-Ras antibody diluted at 1:2000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.



Publications

Ko A, Hasanain M, Oh YT et al. LZTR1 mutation mediates oncogenesis through stabilization of EGFR and AXL
Cancer discovery 2022-11-29 [PMID: 36445254]

Moser R, Gurley KE, Nikolova O et al. Synthetic lethal kinases in Ras/p53 mutant squamous cell carcinoma
Oncogene 2022-05-10 [PMID: 35538224] (WB, Mouse)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-42864

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-50831-0.1mg	Recombinant Human HRAS His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-42864

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

