

Product Datasheet

Myc Epitope Tag Antibody (2) - Azide and BSA Free NBP2-43627

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-43627

Updated 2/23/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-43627




NBP2-43627

Myc Epitope Tag Antibody (2) - Azide and BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2
Preservative	No Preservative
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS
Product Description	
Host	Mouse
Species	Epitope Tag
Immunogen	The immunogen used to generate this antibody corresponds to Myc Epitope Tag
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot 1:5000-1:20000, Immunocytochemistry/ Immunofluorescence 1:100-1:2000, Immunoprecipitation 1:100-1:500

Images


Western Blot: Myc Epitope Tag Antibody (2) [NBP2-43627] - Myc tag antibody [0002] Immunoprecipitation of Myc tagged protein. HeLa transfected lysates were subjected to immunoprecipitation using (B) normal mouse IgG or (C) 2.5 ug of anti-Myc tag antibody [0002] . (A) Input, 20ug of HeLa transfected lysates. The precipitated protein was detected. EasyBlot anti-Mouse IgG Kit was used in Western blot.



Immunocytochemistry/Immunofluorescence: Myc Epitope Tag Antibody (2) [NBP2-43627] - Mock and transfected 293T cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: Myc tag stained by Myc tag antibody [0002] (NBP2-43627) diluted at 1:2000.



Western Blot: Myc Epitope Tag Antibody (2) [NBP2-43627] - Non-transfected (-) and transfected (+) 293T whole cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with Myc Epitope Tag antibody [GT0002] (NBP2-43627) diluted at 1:5000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-43627

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-43627

Earn gift cards/discounts by submitting a publication using this product:

www.novusbio.com/publications

