

Product Datasheet

SERBP1 Antibody (OTI5G2)

NBP2-45763

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-45763

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-45763



NBP2-45763

SERBP1 Antibody (OTI5G2)

Product Information

Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI5G2
Preservative	0.02% Sodium Azide
Isotype	IgG2a
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
Target Molecular Weight	44.1 kDa

Product Description


Host	Mouse
Gene ID	26135
Gene Symbol	SERBP1
Species	Human, Mouse
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions.
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-262 of human SERBP1 (NP_001018078) produced in E.coli.

Product Application Details


Applications	Western Blot, Immunohistochemistry
Recommended Dilutions	Western Blot 1:2000, Immunohistochemistry 1:150

Images


Western Blot: SERBP1 Antibody (5G2) [NBP2-45763] - Analysis of HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SERBP1.




Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] - Analysis of Human endometrium tissue. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120C for 3min)




Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] - Analysis of Carcinoma of Human lung tissue. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120C for 3min)




Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] - Analysis of Human tonsil tissue. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120C for 3min)




Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] - Analysis of Human thyroid tissue. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120C for 3min)




Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] - Analysis of Adenocarcinoma of Human endometrium tissue. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120C for 3min)




Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] - Analysis of Carcinoma of Human prostate tissue. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120C for 3min)




Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] - Analysis of Adenocarcinoma of Human colon tissue. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120C for 3min)




Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] - Analysis of Carcinoma of Human thyroid tissue. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120C for 3min)




Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] -
Analysis of Human pancreas tissue. (Heat-induced epitope retrieval by
10mM citric buffer, pH6.0, 120C for 3min)



Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] -
Analysis of Adenocarcinoma of Human ovary tissue. (Heat-induced
epitope retrieval by 10mM citric buffer, pH6.0, 120C for 3min)




Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] -
Analysis of Human liver tissue. (Heat-induced epitope retrieval by 10mM
citric buffer, pH6.0, 120C for 3min)




Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] -
Analysis of Human prostate tissue. (Heat-induced epitope retrieval by
10mM citric buffer, pH6.0, 120C for 3min)




Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] -
Analysis of Human breast tissue. (Heat-induced epitope retrieval by
10mM citric buffer, pH6.0, 120C for 3min)



Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] -
Analysis of Human Kidney tissue. (Heat-induced epitope retrieval by
10mM citric buffer, pH6.0, 120C for 3min)



Immunohistochemistry: SERBP1 Antibody (5G2) [NBP2-45763] -
Analysis of Carcinoma of Human pancreas tissue. (Heat-induced epitope
retrieval by 10mM citric buffer, pH6.0, 120C for 3min)



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-45763

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)
H00026135-P01-10ug	Recombinant Human SERBP1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-45763

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

