

Product Datasheet

CD1c/BDCA-1 Antibody (OTI2F4) NBP2-46123

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-46123

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-46123



NBP2-46123

CD1c/BDCA-1 Antibody (OTI2F4)


Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI2F4
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
Target Molecular Weight	35.7 kDa

Product Description	
Host	Mouse
Gene ID	911
Gene Symbol	CD1C
Species	Human
Immunogen	Full length human recombinant protein of human CD1C(NP_001756) produced in HEK293T cell.


Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:200 - 1:1000, Immunohistochemistry 1:150, Immunocytochemistry/ Immunofluorescence 1:100, Immunohistochemistry-Paraffin 1:150

Images


Western Blot: CD1c/BDCA-1 Antibody (OTI2F4) [NBP2-46123] - Analysis of HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CD1c .




Immunocytochemistry/Immunofluorescence: CD1c/BDCA-1 Antibody (OTI2F4) [NBP2-46123] - Analysis of COS7 cells transiently transfected by pCMV6-ENTRY CD1C.



Immunohistochemistry: CD1c/BDCA-1 Antibody (OTI2F4) [NBP2-46123] - Analysis of Human endometrium tissue.




Western Blot: CD1c/BDCA-1 Antibody (OTI2F4) [NBP2-46123] - Analysis of Jurkat extracts (10ug).




Immunocytochemistry/Immunofluorescence: CD1c/BDCA-1 Antibody (OTI2F4) [NBP2-46123] - Analysis of HeLa cells.




Immunohistochemistry: CD1c/BDCA-1 Antibody (OTI2F4) [NBP2-46123]
- Analysis of Human lymph node tissue. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120C for 3min)



Immunohistochemistry: CD1c/BDCA-1 Antibody (OTI2F4) [NBP2-46123]
- Analysis of Human tonsil tissue. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120C for 3min)



Immunohistochemistry: CD1c/BDCA-1 Antibody (OTI2F4) [NBP2-46123]
- Analysis of Carcinoma of Human lung tissue.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-46123

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
H00000911-Q01-10ug	Recombinant Human CD1c/BDCA-1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-46123

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

