

Product Datasheet

CD34 Antibody (HPCA1/1171) - Azide and BSA Free NBP2-47911-0.1mg

Unit Size: 0.1 mg

Store at -20 to -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-47911

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-47911



NBP2-47911-0.1mg

CD34 Antibody (HPCA1/1171) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at -20 to -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HPCA1/1171
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	Protein A or G purified
Buffer	10 mM PBS


Product Description	
Description	1.0 mg/ml of antibody purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS WITHOUT BSA & azide. Also available at 200 ug/ml WITH BSA & azide (NBP2-44569). Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80C.
Host	Mouse
Gene ID	947
Gene Symbol	CD34
Species	Human, Rat
Marker	Hematopoietic Stem Cell & Endothelial Marker
Specificity/Sensitivity	This antibody recognizes a carbohydrate epitope on a single chain, transmembrane, heavily glycosylated protein of 90-120kDa, which is identified as CD34 (VI international workshop on human differentiation antigens). Its expression is a hallmark for identifying pluripotent hematopoietic stem or progenitor cells. Its expression is gradually lost as lineage committed progenitors differentiate. CD34 is a marker of choice for staining blasts in acute myeloid leukemia. In addition, it is expressed by soft tissue tumors, such as solitary fibrous tumor and gastrointestinal stromal tumor. CD34 expression is also found in vascular endothelium. Additionally, proliferating endothelial cells overexpress this molecule than the non-proliferating endothelial cells. Anti-CD34 labels 85% of angiosarcoma and Kaposi sarcoma, but shows low specificity.
Immunogen	Recombinant full-length human CD34 protein (Uniprot: P28906)

Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, CyTOF-ready, Immunofluorescence
Recommended Dilutions	Western Blot 0.5 - 1.0 ug/ml, Flow Cytometry 0.5 - 1 ug/million cells in 0.1 ml, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1-2 ug/ml, Immunohistochemistry-Paraffin 2 - 4 ug/ml, Immunofluorescence 0.5 - 1.0 ug/ml, CyTOF-ready
Application Notes	Immunohistochemistry (Formalin-fixed): 2-4ug/ml for 30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95C followed by cooling at RT for 20 minutes. Optimal dilution for a specific application should be determined.




Images


Immunohistochemistry-Paraffin: CD34 Antibody (HPCA1/1171) - Azide and BSA Free [NBP2-47911] - Human Tonsil stained with CD34 Monoclonal Antibody (HPCA1/1171)




Flow Cytometry: CD34 Antibody (HPCA1/1171) - Azide and BSA Free [NBP2-47911] - KG-1 cells using CD34 Monoclonal Antibody (HPCA1/1171) (red) & isotype control (green).



Immunohistochemistry-Paraffin: CD34 Antibody (HPCA1/1171) - Azide and BSA Free [NBP2-47911] - Human Angiosarcoma stained with CD34 Monoclonal Antibody (HPCA1/1171)



Flow Cytometry: CD34 Antibody (HPCA1/1171) - Azide and BSA Free [NBP2-47911] - Percent of EPCs (CD133 + CD34 +) reduced in the XYK group detected by flow cytometry. Data are mean +/- SD analyzed by one-way analysis of variance. **P < 0.01. CSDH, chronic subdural hematoma; EPC, endothelial progenitor cell; XYK, Xiaoyukang Jiaonang. PE conjugated CD34 antibody used (NBP2-47911PE). Image collected and cropped by CiteAb from the following publication (<http://pubmed.ncbi.nlm.nih.gov/32328124/>), licensed under a CC-BY license.



Publications

Wang L, Pan Y, Liu M et al. Wen-Shen-Tong-Luo-Zhi-Tong Decoction regulates bone-fat balance in osteoporosis by adipocyte-derived exosomes *Pharmaceutical biology* 2023-12-01 [PMID: 36999351] (FLOW, Mouse)

Deng W, Jo J, Tanaka T et al. A Senomorphic?Conjugated Scaffold for Application of Senescent Cells in Regenerative Medicine *Advanced Therapeutics* 2023-02-16 (FLOW, Flow)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-47911-0.1mg

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
NBP2-53044-50ug	Recombinant Human CD34 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-47911

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

