

# Product Datasheet

## CD63 Antibody (NKI/C3 + LAMP3/968) [DyLight 488] NBP2-47937G

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-47937G](http://www.novusbio.com/NBP2-47937G)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-47937G](http://www.novusbio.com/reviews/destination/NBP2-47937G)



**NBP2-47937G**

CD63 Antibody (NKI/C3 + LAMP3/968) [DyLight 488]

| <b>Product Information</b>         |  |
|------------------------------------|--|
| <b>Unit Size</b>                   | 0.1 ml   |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.  |
| <b>Storage</b>                     | Store at 4C in the dark.   |
| <b>Clonality</b>                   | Monoclonal   |
| <b>Clone</b>                       | NKI/C3 + LAMP3/968   |
| <b>Preservative</b>                | 0.05% Sodium Azide   |
| <b>Isotype</b>                     | IgG1 Kappa/IgG1 Kappa  |
| <b>Conjugate</b>                   | DyLight 488  |
| <b>Purity</b>                      | Protein A or G purified  |
| <b>Buffer</b>                      | 50mM Sodium Borate   |
| <b>Product Description</b>         |  |
| <b>Host</b>                        | Mouse  |
| <b>Gene ID</b>                     | 967  |
| <b>Gene Symbol</b>                 | CD63   |
| <b>Species</b>                     | Human  |
| <b>Marker</b>                      | Late Endosomes Marker  |
| <b>Specificity/Sensitivity</b>     | This monoclonal antibody recognizes protein of 26kDa-60kDa, which is identified as CD63. The tetraspanins are integral membrane proteins expressed on cell surface and granular membranes of hematopoietic cells and are components of multi-molecular complexes with specific integrins. The tetraspanin CD63 is a lysosomal membrane glycoprotein that translocates to the plasma membrane after platelet activation. CD63 is expressed on activated platelets, monocytes and macrophages, and is weakly expressed on granulocytes, T cell and B cells. It is located on the basophilic granule membranes and on the plasma membranes of lymphocytes and granulocytes. CD63 is a member of the TM4 superfamily of leukocyte glycoproteins that includes CD9, CD37 and CD53, which contain four transmembrane regions. CD63 may play a role in phagocytic and intracellular lysosome-phagosome fusion events. CD63 deficiency is associated with Hermansky-Pudlak syndrome and is strongly expressed during the early stages of melanoma progression. |
| <b>Immunogen</b>                   | Smooth plasma membrane fraction of MeWo cells (NKI/C3); Recombinant human full-length CD63 protein (LAMP3/968)   |
| <b>Product Application Details</b> |  |
| <b>Applications</b>                | Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, CyTOF-ready, Immunofluorescence  |
| <b>Recommended Dilutions</b>       | Flow Cytometry, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, Immunofluorescence, CyTOF-ready  |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.   |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-47937G**

---

|                   |                                      |
|-------------------|--------------------------------------|
| H00000967-G01-2ug | Recombinant Human CD63 Protein       |
| 203-IL-010        | IL-3 [Unconjugated]                  |
| 5048-CD-050       | CD63 [Unconjugated]                  |
| AF3628            | CD31/PECAM-1 Antibody [Unconjugated] |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-47937G](http://www.novusbio.com/reviews/submit/NBP2-47937G)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

