

# Product Datasheet

## Kallikrein 3/PSA Antibody (1A7) [DyLight 594] NBP2-47968DL594

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-47968DL594](http://www.novusbio.com/NBP2-47968DL594)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-47968DL594](http://www.novusbio.com/reviews/destination/NBP2-47968DL594)



**NBP2-47968DL594**

Kallikrein 3/PSA Antibody (1A7) [DyLight 594]

| <b>Product Information</b>         |  |
|------------------------------------|--|
| <b>Unit Size</b>                   | 0.1 ml   |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.  |
| <b>Storage</b>                     | Store at 4C in the dark.   |
| <b>Clonality</b>                   | Monoclonal   |
| <b>Clone</b>                       | 1A7  |
| <b>Preservative</b>                | 0.05% Sodium Azide   |
| <b>Isotype</b>                     | IgG1 Kappa   |
| <b>Conjugate</b>                   | DyLight 594  |
| <b>Purity</b>                      | Protein A or G purified  |
| <b>Buffer</b>                      | 50mM Sodium Borate   |
| <b>Product Description</b>         |  |
| <b>Host</b>                        | Mouse  |
| <b>Gene ID</b>                     | 354  |
| <b>Gene Symbol</b>                 | KLK3   |
| <b>Species</b>                     | Human, Feline  |
| <b>Specificity/Sensitivity</b>     | Recognizes a single protein of 33-34kDa, identified as the prostate specific antigen (PSA). This monoclonal antibody is highly specific to PSA and stains prostatic secretory and ductal epithelium in both normal and neoplastic tissues. PSA is a chymotrypsin-like serine protease (kallikrein family) exclusively produced by the prostate epithelium, and abundant in seminal fluid. PSA can be detected in the sera of patients with prostatic carcinoma. It is predominantly complexed to a liver-derived serine protease inhibitor, alpha-1-antichymotrypsin (ACT). A higher proportion of serum PSA is complexed to ACT in prostate cancer than in benign prostate hyperplasia. This monoclonal antibody makes an excellent pair with monoclonal antibody A67-B/E3 for PSA tests. |
| <b>Immunogen</b>                   | Kallikrein 3/PSA from human sperm plasma (Uniprot: P07288)   |
| <b>Notes</b>                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.  |
| <b>Product Application Details</b> |  |
| <b>Applications</b>                | Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, CyTOF-ready, Immunofluorescence  |
| <b>Recommended Dilutions</b>       | Flow Cytometry, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, Immunofluorescence, CyTOF-ready  |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.   |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-47968DL594**

---

|                 |   |
|-----------------|---|
| NBP1-43319DL594 | Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1) [DyLight 594] |
| NBP2-61380-1mg  | Recombinant Human Kallikrein 3/PSA Protein                  |
| M6000B-1        | IL-6 [HRP]  |
| DKK300          | Kallikrein 3/PSA [HRP]                                      |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-47968DL594](http://www.novusbio.com/reviews/submit/NBP2-47968DL594)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

