

Product Datasheet

TNF-alpha Antibody (52B83) [PE] NBP2-48484PE

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-48484PE

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-48484PE



NBP2-48484PE

TNF-alpha Antibody (52B83) [PE]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	52B83
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	PE
Purity	Protein A or G purified
Buffer	PBS

Product Description	
Description	This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.
Host	Mouse
Gene ID	7124
Gene Symbol	TNF
Species	Human, Mouse, Rabbit, Sheep
Reactivity Notes	Expected cross reactivity based on sequence identity: Guinea pig, Rhesus Monkey, Mouse. Sheep data from customer review.
Specificity/Sensitivity	This detects both natural and recombinant TNFa. It does not cross-react with TNF beta or lymphotoxin. Non-specific background staining is observed in connective tissues.
Immunogen	The monoclonal antibody 52B83 reacts with tumor necrosis factor alpha (TNF-alpha).

Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunohistochemistry, CyTOF-ready
Recommended Dilutions	Western Blot, Flow Cytometry, Immunohistochemistry, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.



Images

TNF-alpha Antibody (52B83) [PE] [NBP2-48484PE] - Vial of PE conjugated antibody. PE has two excitation maxima, 498 nm excited by the Blue laser (488 nm) and 565 nm excited by the Yellow-Green laser (561 nm). Both result in emission at 578 nm.



PE

LASER (nm)	FILTER
Blue (488)	582/15
Y-G (561)	

EXCITATION MAX (nm)	EMISSION MAX (nm)
498	578
565	



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-48484PE

NBP1-97005PE	Mouse IgG1 Isotype Control (MG1) [PE]
NBP2-35076-10ug	Recombinant Human TNF-alpha Protein
M6000B-1	IL-6 [HRP]
210-TA-005	TNF-alpha [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-48484PE

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

