

Product Datasheet

RNA Polymerase II/POLR2A [p Ser5] Antibody (CTD 4H8) [mFluor Violet 500 SE] NBP2-50274MFV500

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-50274MFV500

Updated 7/11/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-50274MFV500



NBP2-50274MFV500

RNA Polymerase II/POLR2A [p Ser5] Antibody (CTD 4H8) [mFluor Violet 500 SE]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	CTD 4H8
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	mFluor Violet 500 SE
Purity	Protein A purified
Buffer	50mM Sodium Borate
Product Description	
Host	Mouse
Gene ID	5430
Gene Symbol	POLR2A
Species	Human, Yeast
Reactivity Notes	Saccharomyces cerevisiae
Specificity/Sensitivity	Please note that this antibody detects both unphosphorylated and phosphorylated forms of the protein.
Immunogen	Peptide: 10 repeats of synthetic peptide YSPTSPS using chemically synthesized phospho-Ser5
Notes	mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.
Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP), Chromatin Immunoprecipitation Sequencing, CyTOF-ready
Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP), CyTOF-ready, Chromatin Immunoprecipitation Sequencing
Application Notes	Optimal dilution of this antibody should be experimentally determined.



Images

RNA Polymerase II/POLR2A [p Ser5] Antibody (CTD 4H8) [mFluor Violet 500 SE] - Vial of mFluor Violet 500 conjugated antibody. mFluor Violet 500 is optimally excited at 410 nm by the Violet laser (405 nm) and has an emission maximum of 501 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP2-50274MFV500

H00005430-Q01-10ug	Recombinant Human RNA Polymerase II/POLR2A GST (N-Term) Protein
H00006908-Q01-10ug	Recombinant Human TATA binding protein TBP GST (N-Term) Protein
NBP1-87785PEP	RNA Polymerase II/POLR2A Recombinant Protein Antigen
NB600-501	beta-Actin Antibody (AC-15)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-50274MFV500

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



