

Product Datasheet

Prohibitin 2 Antibody - BSA Free NBP2-50550

Unit Size: 400 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-50550

Updated 2/23/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-50550



NBP2-50550

Prohibitin 2 Antibody - BSA Free

Product Information	
Unit Size	400 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS

Product Description	
Host	Rabbit
Gene ID	11331
Gene Symbol	PHB2
Species	Human, Mouse
Immunogen	KLH conjugated synthetic peptide between 225-255 amino acids from human PHB2

Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:1000, Immunohistochemistry 1:50 - 1: 100, Immunohistochemistry-Paraffin



Images


Western Blot: Prohibitin 2 Antibody [NBP2-50550] - Analysis using a PHB2 polyclonal antibody in MDA-MB435 cell lysates (35 ug per lane).

95
72
55
36
28




Western Blot: Prohibitin 2 Antibody [NBP2-50550] - Analysis of PHB2 (arrow) using a PHB2 polyclonal antibody in 293 cell lysates (2 ug/lane) either nontransfected (Lane 1) or transiently transfected with the PHB2 gene (Lane 2).

135
95
72
52
42
34
26
(-) (+)



Western Blot: Prohibitin 2 Antibody [NBP2-50550] - Analysis using a PHB2 polyclonal antibody in mouse kidney tissue lysates (35 ug per lane).

95
72
55
36
28





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-50550

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-23365	Recombinant Human Prohibitin 2 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-50550

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

