

Product Datasheet

Recombinant Human pan Cadherin His Protein NBP2-52307-0.05mg

Unit Size: 0.05 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-52307

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-52307



NBP2-52307-0.05mg

Recombinant Human pan Cadherin His Protein

Product Information	
Unit Size	0.05 mg
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Preservative	No Preservative
Purity	>90%, by SDS-PAGE
Buffer	PBS (pH 7.4), 10% glycerol
Target Molecular Weight	76.8 kDa


Product Description	
Description	A recombinant protein with a C-Terminal His-tag and corresponding to the amino acids 24-709 of Human pan Cadherin Source:<i>Baculovirus</i> Amino Acid Sequence: QEPEPCHPGF DAESYTFTVP RRHLERGRVL GRVNFEDCTG RQRTAYFSLD TRFKVGTG DV ITVKRPLRFH NPQIHFLVYA WDSTYRKfst KVTlNTVGHH HRPPPHQASV SGIQAELLTF PNSSPGLRRQ KRDWVIPPIS CPENEKGPFP KNLVQIKSNK DKEGKVFYSI TGQGADTPPV GVFIERETG WLKVTEPLDR ERIATYTLFS HAVSSNGNAV EDPMEILITV TDQNDNKPEF TQEVFKGSVM EGALPGTSVM EVTATDADDD VNTYNAAIAY TILSQDPELP DKNMFTINRN TGVISVTTG LDRESFPTYT LVVQAADLQG EGLSTTATAV ITVTDNDNP PIFNPTTYKG QVPENEANVV ITTLKVTDAD APNTPAWEAV YTLNDDGGQ FVVTNPNVNN DGILKTAKGL DFEAKQQYIL HVAVTNVVPF EVSLTTSTAT VTVDVLDVNE APIFVPPEKR VEVSEDFGVG QEITSYTAQE PDTFMEQKIT YRIWRDTANW LEINPDTGAI STRAELDRED FEHVKNSTYT ALIATDNGS PVATGTGTL LILSDVNDNA PIPEPRTIFF CERNPKPQVI NIIDADLPPN TSPFTAELTH GASANWTIQY NDPTQESIIL KPKMALEVGD YKINLKLMDN QNKDQVTTL VSVCDCEGAA GVCCKAQPVE AGLQIPA
Gene ID	999
Gene Symbol	CDH1
Species	Human
Endotoxin Note	< 1.0 EU per 1 microgram of protein (determined by LAL method)

Product Application Details	
Applications	SDS-Page
Recommended Dilutions	SDS-Page



Images

SDS-Page: Recombinant Human pan Cadherin His Protein [NBP2-52307] - 3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-52307-0.05mg

DVE00	VEGF [HRP]
NB200-592	pan Cadherin Antibody
NB300-223	Vimentin Antibody
NBP1-91258	Fibronectin Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-52307

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

