

Product Datasheet

Gemin 3 Antibody (1C9F6) - BSA Free NBP2-52454

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-52454

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-52454



NBP2-52454

Gemin 3 Antibody (1C9F6) - BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1C9F6
Preservative	0.05% Sodium Azide
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	92.2 kDa

Product Description	
Host	Mouse
Gene ID	11218
Gene Symbol	DDX20
Species	Human, Rat
Immunogen	Purified recombinant fragment of human Gemin 3 (AA: 725-824) expressed in E. Coli.


Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, CyTOF-ready
Recommended Dilutions	Western Blot 1:500-1:2000, Flow Cytometry 1:200-1:400, ELISA 1:10000, Immunocytochemistry/ Immunofluorescence 1:200-1:1000, CyTOF-ready
Application Notes	This antibody is Cytof ready.

Images


Western Blot: Gemin 3 Antibody (1C9F6) [NBP2-52454] - Analysis using DDX20 mAb against human DDX20 (AA: 725-824) recombinant protein. (Expected MW is 38.4 kDa)




Immunocytochemistry/Immunofluorescence: Gemin 3 Antibody (1C9F6) [NBP2-52454] - Analysis of HeLa cells using DDX20 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin.




Flow Cytometry: Gemin 3 Antibody (1C9F6) [NBP2-52454] - Analysis of HeLa cells using DDX20 mouse mAb (green) and negative control (red).




Western Blot: Gemin 3 Antibody (1C9F6) [NBP2-52454] - Western blot analysis using DDX20 mouse mAb against Jurkat (1), Hela (2), NTERA-2 (3), HL7702 (4), K562 (5), and C6 (6) cell lysate.



Western Blot: Gemin 3 Antibody (1C9F6) [NBP2-52454] - Analysis using DDX20 mAb against HEK293 (1) and DDX20 (AA: 725-824)-hIgGfc transfected HEK293 (2) cell lysate.



ELISA: Gemin 3 Antibody (1C9F6) [NBP2-52454] - Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-52454

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
H00011218-P02-10ug	Recombinant Human Gemin 3 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-52454

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

