

Product Datasheet

His Tag Antibody (3D5) - Chimeric - Azide and BSA Free NBP2-52637-0.2mg

Unit Size: 0.2 mg

Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.

www.novusbio.com



technical@novusbio.com

Reviews: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-52637

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-52637



NBP2-52637-0.2mg

His Tag Antibody (3D5) - Chimeric - Azide and BSA Free

Product Information	
Unit Size	0.2 mg
Concentration	1 mg/ml
Storage	Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.
Clonality	Monoclonal
Clone	3D5
Preservative	0.02% Proclin 300
Isotype	IgG Kappa
Purity	Protein A purified
Buffer	PBS


Product Description	
Host	Rabbit
Species	Epitope Tag
Specificity/Sensitivity	This antibody binds to carboxy-terminal oligo-histidine tags (6x His tags) of fusion proteins. It also has affinity for internal His-tags but no affinity for N-terminal His-tags.
Immunogen	This His Tag Antibody (3D5) - Chimeric was developed against purified his tagged Fab fragment.

Product Application Details	
Applications	Western Blot, Simple Western, ELISA
Recommended Dilutions	Western Blot 1:100 - 1:2000, Simple Western, ELISA 1:100 - 1:2000
Application Notes	This is a chimeric rabbit antibody, based on the original mouse IgG2b format, created for improved compatibility with existing reagents, assays and techniques.




Images

Western Blot: His Tag Antibody (3D5) [NBP2-52637] - Staining of Multitag protein at 25 ng (A), 10 ng (B) and 5 ng (C) amounts. Used at a dilution of 1 ug/mL. Primary incubation was 1 hour. Detected by chemiluminescence.



Simple Western: His Tag Antibody (3D5) - Chimeric [NBP2-52637] - His Tag Antibody (3D5) [NBP2-52637] - Human purified protein. Simple Western image submitted by a verified customer review.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-52637-0.2mg

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-52637

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

