

Product Datasheet

Fc gamma RII/RIII (CD32/CD16) Antibody (2.4G2) **[PE/Cy5.5]** **NBP2-52644PECY55**

Unit Size: 0.1 ml

Store at 4C in the dark. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-52644PECY55

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-52644PECY55



NBP2-52644PECY55

Fc gamma RII/RIII (CD32/CD16) Antibody (2.4G2) [PE/Cy5.5]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark. Do not freeze.
Clonality	Monoclonal
Clone	2.4G2
Preservative	0.05% Sodium Azide
Isotype	IgG Kappa
Conjugate	PE/Cy5.5
Purity	Protein A purified
Buffer	PBS

Product Description	
Host	Rabbit
Gene ID	2214
Gene Symbol	FCGR3A
Species	Mouse
Reactivity Notes	Mouse
Specificity/Sensitivity	This antibody binds to an epitope on the extracellular domains of the mouse Fc (gamma)III and Fc(gamma)II, thereby inhibiting their ability to bind antibodies via their Fc region. receptors. Also binding the Fc(gamma)I receptor (CD64) via its Fc domain.
Immunogen	Mouse BALB/c Macrophage J774 Cell Line.

Product Application Details	
Applications	Flow Cytometry, Product Image
Recommended Dilutions	Flow Cytometry, Product Image
Application Notes	Optimal dilution of this antibody should be experimentally determined. For optimal results using our Tandem dyes, please avoid prolonged exposure to light or extreme temperature fluctuations. These can lead to irreversible degradation or decoupling. When staining intracellular targets, specific attention to the fixation and permeabilization steps in your flow protocol may be required. Please contact our technical support team at technical@novusbio.com if you have any questions.



Images

Product Image: CD16/CD32 Antibody (2.4G2) [PE/Cy5.5] [NBP2-52644PECY55] - Vial of PE/Cy5.5 conjugated antibody. PE/Cy5.5 has two excitation maxima, 498 nm excited by the Blue laser (488 nm) and 566 nm excited by the Yellow-Green laser (561 nm). Both result in emission at 700 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-52644PECY55

210-TA-005	TNF-alpha [Unconjugated]
202-IL-010	IL-2 [Unconjugated]
DC140	CD14 [HRP]
NB110-89474	CD11b Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-52644PECY55

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

