

Product Datasheet

Vimentin Antibody (CL0157) - BSA Free NBP2-52867

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-52867

Updated 2/27/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-52867



NBP2-52867

Vimentin Antibody (CL0157) - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL0157
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH 7.2) and 40% Glycerol
Target Molecular Weight	53.6 kDa

Product Description	
Host	Mouse
Gene ID	7431
Gene Symbol	VIM
Species	Human
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions
Immunogen	This antibody was developed using a recombinant protein derived from P08670, with the exact immunogen sequence remaining proprietary.


Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin, Knockout Validated
Recommended Dilutions	Western Blot 1 ug/ml, Immunohistochemistry 1:2500 - 1:5000, Immunohistochemistry-Paraffin 1:2500 - 1:5000, Knockout Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images


Western Blot: Vimentin Antibody (CL0157) [NBP2-52867] - Western blot shows lysates of K562 human Chronic Myelogenous Leukemia parental cell line and Vimentin knockout (KO) K562 cell line. PVDF membrane was probed with 0.7 ug/mL of Mouse Anti-Human Vimentin Monoclonal Antibody (Catalog # NBP2-52867) followed by HRP-conjugated Anti-Mouse IgG Secondary Antibody (Catalog #HAF018). Specific band was detected for Vimentin at approximately 55 kDa (as indicated) in the parental K562 cell line, but is not detectable in the knockout K562 cell line. This experiment was conducted under reducing conditions.



Western Blot: Vimentin Antibody (CL0157) [NBP2-52867] - Analysis in human cell lines U-251MG and MCF-7 using Anti-VIM antibody. Corresponding VIM RNA-seq data are presented for the same cell lines. Loading control: Anti-GAPDH.




Immunohistochemistry-Paraffin: Vimentin Antibody (CL0157) [NBP2-52867] - Staining of human skin shows distinct positivity in fibroblasts.



Western Blot: Vimentin Antibody (CL0157) [NBP2-52867] - Lane 1: Marker [kDa] Lane 2: Human tonsil tissue lysate




Western Blot: Vimentin Antibody (CL0157) [NBP2-52867] - Lane 1: Marker [kDa], Lane 2: Human cell line HeLa cytoplasmic fraction, Lane 3: Human cell line HeLa membrane fraction, Lane 4: Human cell line HeLa nuclear fraction, Lane 5: Human cell line HeLa chromatin fraction, Lane 6: Human cell line HeLa cytosol




Immunohistochemistry-Paraffin: Vimentin Antibody (CL0157) [NBP2-52867] - Staining of human cerebellum shows immunoreactivity in the endothelial cells.




Immunohistochemistry-Paraffin: Vimentin Antibody (CL0157) [NBP2-52867] - Staining of human kidney shows strong immunoreactivity in renal glomeruli but not tubuli cells.



Immunohistochemistry-Paraffin: Vimentin Antibody (CL0157) [NBP2-52867] - Staining of human prostate shows strong immunoreactivity in the connective tissue.



Immunohistochemistry-Paraffin: Vimentin Antibody (CL0157) [NBP2-52867] - Staining of human rectum shows strong positivity in mesenchymal and lymphoid cells.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-52867

NBP2-52867PEP	Vimentin Recombinant Protein Antigen
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-52867

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

