

Product Datasheet

TAZ/WWTR1 Antibody (CL0371) - BSA Free NBP2-52926

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-52926

Updated 2/27/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-52926



NBP2-52926

TAZ/WWTR1 Antibody (CL0371) - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL0371
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Host	Mouse
Gene ID	25937
Gene Symbol	WWTR1
Species	Human
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions
Immunogen	This antibody was developed using a recombinant protein derived from Q9GZV5, with the exact immunogen sequence remaining proprietary.

Product Application Details


Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Knockdown Validated
Recommended Dilutions	Western Blot 1 ug/ml, Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 2-10 ug/ml, Immunohistochemistry-Paraffin 1:50 - 1:200, Knockdown Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.

Images


Immunohistochemistry-Paraffin: TAZ/WWTR1 Antibody (CL0371) [NBP2-52926] - Staining in human placenta and lymph node tissues. Corresponding WWTR1 RNA-seq data are presented for the same tissues.




Western Blot: TAZ/WWTR1 Antibody (CL0371) [NBP2-52926] - Analysis in EFO-21 cells transfected with control siRNA, target specific siRNA probe #1 and #2, using Anti-WWTR1 antibody. Remaining relative intensity is presented. Loading control: Anti-PPIB.




Immunocytochemistry/Immunofluorescence: TAZ/WWTR1 Antibody (CL0371) [NBP2-52926] - Staining in SKMEL cell line with Anti-WWTR1 monoclonal antibody, showing spotty nuclear staining in green. Microtubule probes are visualized in red (where available). Antibody staining is shown in green.




Western Blot: TAZ/WWTR1 Antibody (CL0371) [NBP2-52926] - Lane 1: Marker [kDa], Lane 2: Human cell line U-251 MG




Immunohistochemistry-Paraffin: TAZ/WWTR1 Antibody (CL0371) [NBP2-52926] - Staining of human breast cancer shows moderate to strong nuclear positivity in tumor cells.




Immunohistochemistry-Paraffin: TAZ/WWTR1 Antibody (CL0371) [NBP2-52926] - Staining of human placenta shows strong nuclear and some cytoplasmic positivity in endothelial cells.



Immunohistochemistry-Paraffin: TAZ/WWTR1 Antibody (CL0371) [NBP2-52926] - Staining of human kidney shows strong nuclear and moderate cytoplasmic positivity in cells in tubules and in glomeruli.



Immunohistochemistry-Paraffin: TAZ/WWTR1 Antibody (CL0371) [NBP2-52926] - Staining of human lymph node shows no nuclear positivity in lymphoid cells as expected.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-52926

NBP2-52926PEP	TAZ/WWTR1 Recombinant Protein Antigen
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-52926

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

