

Product Datasheet

Histone Deacetylase 4/HDAC4 [p Ser632] Antibody - BSA Free NBP2-54636-0.05mg

Unit Size: 0.05 mg

Store at -20 degrees C. Avoid freeze/thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-54636

Updated 10/23/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-54636



NBP2-54636-0.05mg


Histone Deacetylase 4/HDAC4 [p Ser632] Antibody - BSA Free

Product Information	
Unit Size	0.05 mg
Concentration	1 mg/ml
Storage	Store at -20 degrees C. Avoid freeze/thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Epitope affinity purified
Buffer	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, and 50% glycerol.
Target Molecular Weight	140 kDa
Product Description	
Host	Rabbit
Gene ID	9759
Gene Symbol	HDAC4
Species	Human, Mouse, Rat
Specificity/Sensitivity	This antibody detects endogenous level of total Histone Deacetylase 4/HDAC4 protein.
Immunogen	HDAC4 (Ab-632) antibody was raised against a peptide sequence around aa. 630~634 (A-Q-S-S-P) derived from Human HDAC4.
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:500-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:200




Images

Western Blot: Histone Deacetylase 4/HDAC4 [p Tyr410] Antibody [NBP2-54636] - Western blot analysis of lysed extracts from HepG2 cells using 4/HDAC4 [p Tyr410].



Immunocytochemistry/Immunofluorescence: Histone Deacetylase 4/HDAC4 [p Tyr410] Antibody [NBP2-54636] - Immunofluorescence staining of methanol-fixed HeLa cells using 4/HDAC4 [p Tyr410]



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-54636-0.05mg

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00009759-P01-10ug	Recombinant Human Histone Deacetylase 4/HDAC4 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-54636

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

