

Product Datasheet

Gelsolin/GSN Antibody - BSA Free NBP2-55565

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-55565

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-55565



NBP2-55565

Gelsolin/GSN Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Description	Novus Biologicals Rabbit Gelsolin/GSN Antibody - BSA Free (NBP2-55565) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	2934
Gene Symbol	GSN
Species	Human
Immunogen	This antibody was developed against a recombinant protein corresponding to the following amino acid sequence: AFVLKTPSAAYLWVG TGASEAEKTGAQELLRLRAQPVQVAEGSEPDGFWEA LGGKAA YRTSPRLKDKKMDAHPPR

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:50 - 1:200, Immunohistochemistry-Paraffin 1:50 - 1:200
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images


Western Blot: Gelsolin/GSN Antibody [NBP2-55565] - Analysis in human cell line RT-4, human cell line U-251 MG and human plasma.




Immunohistochemistry-Paraffin: Gelsolin/GSN Antibody [NBP2-55565] - Staining of human liver shows moderate to strong cytoplasmic positivity in Kupffer cells.




Immunohistochemistry-Paraffin: Gelsolin/GSN Antibody [NBP2-55565] - Immunohistochemical staining of human tonsil shows strong cytoplasmic positivity in dendritic cells.



Immunohistochemistry-Paraffin: Gelsolin/GSN Antibody [NBP2-55565] - Staining of human cervix, uterine shows moderate to strong cytoplasmic positivity in Langerhans cells.



Immunohistochemistry-Paraffin: Gelsolin/GSN Antibody [NBP2-55565] - Staining of human kidney shows moderate to strong cytoplasmic positivity in cells in tubules.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-55565

NBP2-55565PEP	Gelsolin/GSN Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-55565

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

