

Product Datasheet

UPF3B Antibody - BSA Free

NBP2-55617

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-55617

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-55617



NBP2-55617

UPF3B Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol


Product Description	
Description	Novus Biologicals Rabbit UPF3B Antibody - BSA Free (NBP2-55617) is a polyclonal antibody validated for use in WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	65109
Gene Symbol	UPF3B
Species	Human
Reactivity Notes	Mouse 80%, Rat 80%
Immunogen	This antibody was developed against a recombinant protein corresponding to the following amino acid sequence: AKKLDKENLSDERASGQSCTLPKRSDSELKDEKPKRPEDESGRDYREREREY ERDQERILRERERLKRQEEERRRQKERYEKEKTFKRKEEEMKKEKDTLRDKG KKAESTESIGSSEKTEKKEEVVKRDRI RNKDR

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml
Application Notes	Immunocytochemistry/Immunofluorescence Fixation Permeabilization: Use PFA/Triton X-100.




Images


Western Blot: UPF3B Antibody [NBP2-55617] - Analysis using Anti-UPF3B antibody NBP2-55617 (A) shows similar pattern to independent antibody NBP1-83134 (B).




Immunocytochemistry/Immunofluorescence: UPF3B Antibody [NBP2-55617] - Staining of human cell line SiHa shows localization to nucleus, nucleoli & cytosol.



Western Blot: UPF3B Antibody [NBP2-55617] - Western blot analysis in human cell line RT-4.



Western Blot: UPF3B Antibody [NBP2-55617] - Analysis in human cell line NTERA-2.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-55617

NBP2-55617PEP	UPF3B Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-55617

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

