

# Product Datasheet

## IGFBP-2 Recombinant Protein Antigen NBP2-57978PEP

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-57978PEP](http://www.novusbio.com/NBP2-57978PEP)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-57978PEP](http://www.novusbio.com/reviews/destination/NBP2-57978PEP)



**NBP2-57978PEP**

## IGFBP-2 Recombinant Protein Antigen

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Preservative	No Preservative
Purity	>80% by SDS-PAGE and Coomassie blue staining
Buffer	PBS and 1M Urea, pH 7.4.
Target Molecular Weight	24 kDa
Product Description	
Description	<p>A recombinant protein antigen with a N-terminal His6-ABP tag corresponding to human IGFBP-2.</p> <p><b>Source:</b> <i>E. coli</i></p> <p><b>Amino Acid Sequence:</b> GLRCYPHPGSELPLQALVMGEGTCEKRRDAEYGASPEQVADNGDDHSEGGL VENHVDS</p> <p><b>Fusion Tag:</b> N-terminal His6ABP (ABP = Albumin Binding Protein derived from Streptococcal Protein G)</p> <p><b>This product is intended to be used as a blocking antigen for antibody competition assays. Any other use of this antigen is done at the risk of the user. The use of this product for commercial production is strictly prohibited. Please contact technical support if you have any questions.</b></p>
Gene ID	3485
Gene Symbol	IGFBP2
Species	Human
Product Application Details	
Applications	Antibody Competition
Recommended Dilutions	Antibody Competition 10 - 100 molar excess
Application Notes	<p>This recombinant antigen is only intended to be used as a blocking agent to confirm antibody specificity with the corresponding antibody, catalog number NBP2-57978.</p> <p>It is purified by IMAC chromatography, and the expected concentration is greater than 0.5 mg/ml.</p> <p>For current lot information, including availability, please contact our technical support team click <a href="mailto:nb-technical@bio-techne.com">nb-technical@bio-techne.com</a></p> <p>For further blocking tide related information and a protocol, click <a href="#">here</a>.</p>



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-57978PEP**

---

236-EG-200	EGF [Unconjugated]
AF674	IGFBP-2 Antibody [Unconjugated]
M6000B-1	IL-6 [HRP]
291-G1-200	IGF-I/IGF-1 [Unconjugated]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-57978PEP](http://www.novusbio.com/reviews/submit/NBP2-57978PEP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

