

Product Datasheet

IKIP Antibody - BSA Free NBP2-58695

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-58695

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-58695



NBP2-58695

IKIP Antibody - BSA Free


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Description	Novus Biologicals Rabbit IKIP Antibody - BSA Free (NBP2-58695) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	121457
Gene Symbol	IKBIP
Species	Human
Reactivity Notes	Mouse 85%, Rat 84%
Immunogen	This antibody was developed against a recombinant protein corresponding to the following amino acid sequence: DSTLRNIRTVKRQEEEDLLRVEEQLGSDTKAIEKLEEEQHALFARDEDLTNKLS DYEPKVEECKTHLPTIESAI


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:1000 - 1:2500, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:1000 - 1:2500
Application Notes	ICC/IF Fixation Permeabilization: Use PFA/Triton X-100. IHC-P Retrieval method: HIER pH6.

Images


Immunohistochemistry-Paraffin: IKIP Antibody [NBP2-58695] - Staining in human placenta and skeletal muscle tissues using anti-IKBIP antibody. Corresponding IKBIP RNA-seq data are presented for the same tissues.




Immunohistochemistry-Paraffin: IKIP Antibody [NBP2-58695] - Staining of human colon, liver, placenta and skeletal muscle using Anti-IKBIP antibody NBP2-58695 (A) shows similar protein distribution across tissues to independent antibody NBP1-92022 (B).




Western Blot: IKIP Antibody [NBP2-58695] - Analysis using Anti-IKBIP antibody NBP2-58695 (A) shows similar pattern to independent antibody NBP1-92022 (B).




Immunocytochemistry/Immunofluorescence: IKIP Antibody [NBP2-58695] - Staining of human cell line U-2 OS shows localization to endoplasmic reticulum. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: IKIP Antibody [NBP2-58695] - Staining of human liver.




Immunohistochemistry-Paraffin: IKIP Antibody [NBP2-58695] - Staining of human placenta shows high expression.



Immunohistochemistry-Paraffin: IKIP Antibody [NBP2-58695] - Staining of human skeletal muscle shows low expression as expected.



Immunohistochemistry-Paraffin: IKIP Antibody [NBP2-58695] - Staining of human colon.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-58695

NBP2-58695PEP	IKIP Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-58695

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

