

# Product Datasheet

## Histone H3 [ac Lys36] Antibody - BSA Free NBP2-59273

Unit Size: 50 ug

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-59273](http://www.novusbio.com/NBP2-59273)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-59273](http://www.novusbio.com/reviews/destination/NBP2-59273)



**NBP2-59273**


Histone H3 [ac Lys36] Antibody - BSA Free

| Product Information         |                                                                                                                                                                                                                       |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 50 ug                                                                                                                                                                                                                 |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                           |
| Storage                     | Store at -20C. Avoid freeze-thaw cycles.                                                                                                                                                                              |
| Clonality                   | Polyclonal                                                                                                                                                                                                            |
| Preservative                | 0.05% Sodium Azide and 0.05% ProClin 300                                                                                                                                                                              |
| Isotype                     | IgG                                                                                                                                                                                                                   |
| Purity                      | Peptide affinity purified                                                                                                                                                                                             |
| Buffer                      | PBS                                                                                                                                                                                                                   |
| Target Molecular Weight     | 15 kDa                                                                                                                                                                                                                |
| Product Description         |                                                                                                                                                                                                                       |
| Description                 | Novus Biologicals Rabbit Histone H3 [ac Lys36] Antibody - BSA Free (NBP2-59273) is a polyclonal antibody validated for use in WB, ELISA and ChIP. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host                        | Rabbit                                                                                                                                                                                                                |
| Gene ID                     | 126961                                                                                                                                                                                                                |
| Gene Symbol                 | H3C14                                                                                                                                                                                                                 |
| Species                     | Human                                                                                                                                                                                                                 |
| Immunogen                   | The exact sequence of the immunogen to this Histone H3 [ac Lys36] antibody is proprietary.                                                                                                                            |
| Product Application Details |                                                                                                                                                                                                                       |
| Applications                | Western Blot, Dot Blot, ELISA, Chromatin Immunoprecipitation (ChIP), Chromatin Immunoprecipitation Sequencing                                                                                                         |
| Recommended Dilutions       | Western Blot 1:1000, ELISA 1:1000, Dot Blot 1:10000, Chromatin Immunoprecipitation (ChIP) 2 ug, Chromatin Immunoprecipitation Sequencing                                                                              |




## Images


**Western Blot: Histone H3 [ac Lys36] Antibody [NBP2-59273]** - Western blot was performed on whole cell (25 ug, lane 1) and histone extracts (15 ug, lane 2) from HeLa cells, and on 1 ug of recombinant histone H2A, H2B, H3 and H4 (lane 3, 4, 5 and 6, respectively) using the antibody against H3K36ac. The antibody was diluted 1:1,000 in TBS- Tween containing 5% skimmed milk. Observed molecular weight is ~16 kDa.




**Chromatin Immunoprecipitation: Histone H3 [ac Lys36] Antibody [NBP2-59273]** - ChIP assays were performed using human HeLa cells, the antibody against H3K36ac and optimized PCR primer sets for qPCR. ChIP was performed using sheared chromatin from 1.5 million cells. A titration of the antibody consisting of 1, 2, 5 and 10 ug per ChIP experiment was analysed. IgG (2 ug/IP) was used as negative IP control. QPCR was performed with primers for a region approximately 1 kb upstream of the ACTB promoter and for the GAPDH promoter, used as positive controls, and for the coding region of the inactive MYOD1 gene and the Sat2 satellite repeat, used as negative controls. Figure shows the recovery, expressed as a % of input (the relative amount of immunoprecipitated DNA compared to input DNA after qPCR analysis).



**Dot Blot: Histone H3 [ac Lys36] Antibody [NBP2-59273]** - To test the cross reactivity of the antibody against H3K36ac, a Dot Blot analysis was performed with peptides containing other histone modifications and the unmodified H3K36. One hundred to 0.2 pmol of the respective peptides were spotted on a membrane. The antibody was used at a dilution of 1:10,000. Figure shows a high specificity of the antibody for the modification of interest.



**ELISA: Histone H3 [ac Lys36] Antibody [NBP2-59273]** - To determine the titer of the antibody, an ELISA was performed using a serial dilution of the antibody directed against H3K36ac in antigen coated wells. The antigen used was a peptide containing the histone modification of interest. By plotting the absorbance against the antibody dilution, the titer of the antibody was estimated to be 1:10,000.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-59273**

---

|             |                                                     |
|-------------|-----------------------------------------------------|
| NBP2-33376H | Blue Marker Antibody (6F4-F6) [HRP]                 |
| HAF008      | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160      | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891  | Rabbit IgG Isotype Control                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-59273](http://www.novusbio.com/reviews/submit/NBP2-59273)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

