

# Product Datasheet

## SCN3B Antibody (S396-29) [Alexa Fluor® 700] NBP2-59318AF700

Unit Size: 100 ul

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-59318AF700](http://www.novusbio.com/NBP2-59318AF700)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-59318AF700](http://www.novusbio.com/reviews/destination/NBP2-59318AF700)



**NBP2-59318AF700**

SCN3B Antibody (S396-29) [Alexa Fluor® 700]

| Product Information         |                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C in the dark.                                                                    |
| Clonality                   | Monoclonal                                                                                  |
| Clone                       | S396-29                                                                                     |
| Preservative                | 0.05% Sodium Azide                                                                          |
| Isotype                     | IgG2b                                                                                       |
| Conjugate                   | Alexa Fluor 700                                                                             |
| Purity                      | Protein G purified                                                                          |
| Buffer                      | 50mM Sodium Borate                                                                          |
| Product Description         |                                                                                             |
| Host                        | Mouse                                                                                       |
| Gene ID                     | 55800                                                                                       |
| Gene Symbol                 | SCN3B                                                                                       |
| Species                     | Human, Mouse, Rat                                                                           |
| Specificity/Sensitivity     | Detects 40kDa. Does not cross-react with NavBeta1, Navbeta2, or Navbeta4.                   |
| Immunogen                   | Fusion protein amino acids 1-215 (full-length) of rat NavBeta3                              |
| Product Application Details |                                                                                             |
| Applications                | Simple Western, Immunocytochemistry/ Immunofluorescence                                     |
| Recommended Dilutions       | Simple Western, Immunocytochemistry/ Immunofluorescence                                     |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-59318AF700**

---

|                    |                                                         |
|--------------------|---------------------------------------------------------|
| NBP2-27231AF700    | Mouse IgG2b Isotype Control (MPC-11) [Alexa Fluor® 700] |
| NBP2-59531-50ug    | Recombinant Human SCN3B His Protein                     |
| H00006327-P01-10ug | Recombinant Human Scn2b GST (N-Term) Protein            |
| NBL1-15743         | SCN3B Overexpression Lysate                             |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-59318AF700](http://www.novusbio.com/reviews/submit/NBP2-59318AF700)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

