

Product Datasheet

Peroxiredoxin 3 Antibody (1F8) - BSA Free **NBP2-59425-100ul**

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-59425

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-59425



NBP2-59425-100ul

Peroxiredoxin 3 Antibody (1F8) - BSA Free

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1F8
Preservative	0.02% Sodium Azide
Isotype	IgG3 Kappa
Purity	Protein A purified
Buffer	PBS (pH 7.4), 10% Glycerol


Product Description	
Description	Novus Biologicals Mouse Peroxiredoxin 3 Antibody (1F8) - BSA Free (NBP2-59425) is a monoclonal antibody validated for use in WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	10935
Gene Symbol	PRDX3
Species	Human
Immunogen	Recombinant human Peroxiredoxin3 (63-256aa) purified from E. coli

Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:1000, ELISA, Immunocytochemistry/ Immunofluorescence 1:100




Images

Western Blot: Peroxiredoxin 3 Antibody (1F8) [NBP2-59425] - Lane 1: Recombinant protein, Lane 2: HeLa cell lysates, Lane 3: HepG2 cell lysates, Lane 4: TF1 cell lysates, Lane 5: U87MG cell lysates, Lane 6: Raji cell lysates, Lane 7: 293T cell lysates, Lane 8: Jurkat cell lysates, Lane 9: MCF7 cell lysates



Immunocytochemistry/Immunofluorescence: Peroxiredoxin 3 Antibody (1F8) [NBP2-59425] - Analysis of Peroxiredoxin3 in HeLa cells line, stained with DAPI (Blue) for nucleus staining and monoclonal anti-human Peroxiredoxin3 antibody (1:100) with goat anti-mouse IgG-Alexa fluor 488 conjugate (Green).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-59425-100ul

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96978	Mouse IgG3 Kappa Light Chain Isotype Control (MG3K)
NBC1-27993	Recombinant Human Peroxiredoxin 3 Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-59425

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

