

Product Datasheet

Peroxiredoxin 6 Antibody (22E7) - BSA Free NBP2-59436-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-59436

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-59436



NBP2-59436-100ul


Peroxiredoxin 6 Antibody (22E7) - BSA Free

| Product Information | |
|-----------------------------|---|
| Unit Size | 100 ul |
| Concentration | 1 mg/ml |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | 22E7 |
| Preservative | 0.02% Sodium Azide |
| Isotype | IgG2a Kappa |
| Purity | Protein A purified |
| Buffer | PBS (pH 7.4), 10% Glycerol |
| Product Description | |
| Host | Mouse |
| Gene ID | 9588 |
| Gene Symbol | PRDX6 |
| Species | Human |
| Immunogen | Recombinant human Peroxiredoxin 6 (1-224aa) purified from E. coli |
| Product Application Details | |
| Applications | Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence |
| Recommended Dilutions | Western Blot 1:1000, Flow Cytometry, ELISA, Immunocytochemistry/ Immunofluorescence 1:1000 |




Images

Western Blot: Peroxiredoxin 6 Antibody (22E7) [NBP2-59436] - The recombinant human PRDX6 protein (50ng), Cell and Mouse tissue lysates (40ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human Peroxiredoxin 6 antibody (1:1000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system. Lane 1.: Recombinant PRDX6 protein Lane 2.: HeLa cell lysate Lane 3.: 293T cell lysate Lane 4.: HepG2 cell lysate Lane 5.: U87MG cell lysate Lane 6.: Ramos cell lysate



Immunocytochemistry/Immunofluorescence: Peroxiredoxin 6 Antibody (22E7) [NBP2-59436] - Analysis of PRDX6 in HeLa cells line, stained with DAPI (Blue) for nucleus staining and monoclonal anti-human PRDX6 antibody (1:100) with goat anti-mouse IgG-Alexa fluor 488 conjugate (Green).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-59436-100ul

| | |
|------------------|---|
| HAF007 | Goat anti-Mouse IgG Secondary Antibody [HRP] |
| NB720-B | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-96981-0.5mg | Mouse IgG2a Kappa Isotype Control (M2AK) |
| NBP1-30188 | Recombinant Human Peroxiredoxin 6 His Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-59436

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

