

# Product Datasheet

## Human TDP1 ELISA Kit (Colorimetric) NBP2-60608

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-60608](http://www.novusbio.com/NBP2-60608)

Updated 5/4/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-60608](http://www.novusbio.com/reviews/destination/NBP2-60608)



**NBP2-60608****Human TDP1 ELISA Kit (Colorimetric)****Product Information**

|                      |                                                                                                          |
|----------------------|----------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 1 Kit                                                                                                    |
| <b>Concentration</b> | Concentration is not relevant for this product. Please see the protocols for proper use of this product. |
| <b>Storage</b>       | Storage of components varies. See protocol for specific instructions.                                    |

**Product Description**


|                             |                                                                                                                |
|-----------------------------|----------------------------------------------------------------------------------------------------------------|
| <b>Description</b>          | This assay employs a quantitative enzyme immunoassay technique that measures the specified antigen in samples. |
| <b>Gene ID</b>              | 55775                                                                                                          |
| <b>Gene Symbol</b>          | TDP1                                                                                                           |
| <b>Species</b>              | Human                                                                                                          |
| <b>Standard Curve Range</b> | 1.563 - 100 ng/mL (example only; lot dependent)                                                                |
| <b>Sensitivity</b>          | 1.3 ng/mL (example only; lot dependent)                                                                        |
| <b>Inter Assay</b>          | %CV < 11.6% (example only; lot dependent)                                                                      |
| <b>Intra Assay</b>          | %CV < 5.3% (example only; lot dependent)                                                                       |
| <b>Assay Type</b>           | Sandwich ELISA                                                                                                 |
| <b>Spiking Recovery</b>     | 87 - 110% (example only; lot dependent)                                                                        |
| <b>Suitable Sample Type</b> | Serum, Plasma, Cell Culture                                                                                    |
| <b>Sample Volume</b>        | 50 uL                                                                                                          |

**Product Application Details**

|                              |                       |
|------------------------------|-----------------------|
| <b>Applications</b>          | ELISA, Sandwich ELISA |
| <b>Recommended Dilutions</b> | ELISA, Sandwich ELISA |

**Images**

ELISA: Human TDP1 ELISA Kit (Colorimetric) [NBP2-60608] - Standard curve.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-60608**

---

|                 |                                         |
|-----------------|-----------------------------------------|
| NBP2-59661-10ug | Recombinant Human TDP1 His Protein      |
| NB100-116       | APE Antibody (13B8E5C2) - BSA Free      |
| NB100-81642     | TDP1 Antibody                           |
| NB100-1607      | Myelin Protein Zero Antibody - BSA Free |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-60608](http://www.novusbio.com/reviews/submit/NBP2-60608)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

