

# Product Datasheet

## CD8 beta Antibody (PG164A)

### NBP2-60955

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-60955](http://www.novusbio.com/NBP2-60955)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-60955](http://www.novusbio.com/reviews/destination/NBP2-60955)



**NBP2-60955**

CD8 beta Antibody (PG164A)

**Product Information**

|                      |                                                                                        |
|----------------------|----------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 mg                                                                                 |
| <b>Concentration</b> | This product is unpurified. The exact concentration of antibody is not quantifiable.   |
| <b>Storage</b>       | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| <b>Clonality</b>     | Monoclonal                                                                             |
| <b>Clone</b>         | PG164A                                                                                 |
| <b>Preservative</b>  | 0.09% Sodium Azide                                                                     |
| <b>Isotype</b>       | IgG2a                                                                                  |
| <b>Purity</b>        | Unpurified                                                                             |
| <b>Buffer</b>        | Ascites                                                                                |

**Product Description**


|                    |                                                                                                                                                                                                                                                                                                                                                                                                 |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | This monoclonal antibody is made in ascites form using mice. The ascites is collected and centrifuged to remove cells and large particles. The ascitic fluid is purified further by filtration through sequential filters with the last filter at 2 um. The concentration of the monoclonal antibody present in the ascites is determined with commercial radial immunodiffusion kits (1mg/ml). |
| <b>Host</b>        | Mouse                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Gene ID</b>     | 926                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Gene Symbol</b> | CD8B                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Species</b>     | Porcine                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen</b>   | Pig mononuclear cells                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Notes</b>       | This monoclonal antibody is made in ascites form using mice. The ascites is collected and centrifuged to remove cells and large particles. The ascitic fluid is purified further by filtration through sequential filters with the last filter at 2 um. The concentration of the monoclonal antibody present in the ascites is determined with commercial radial immunodiffusion kits (1mg/ml). |

**Product Application Details**

|                              |                              |
|------------------------------|------------------------------|
| <b>Applications</b>          | Flow Cytometry               |
| <b>Recommended Dilutions</b> | Flow Cytometry 1:10 - 1:1000 |

**Images**

Flow Cytometry: CD8 beta Antibody (PG164A) [NBP2-60955] - Obtained with pig leukocytes





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

[nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

[canada.inquires@bio-techne.com](mailto:canada.inquires@bio-techne.com)

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

[info.EMEA@bio-techne.com](mailto:info.EMEA@bio-techne.com)

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)

Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)

Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)

General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-60955**

|                  |                                                         |
|------------------|---------------------------------------------------------|
| HAF007           | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B          | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-96778       | Mouse IgG2a Isotype Control (M2A)                       |
| NBP2-51536-0.1mg | Recombinant Human CD8 beta His Protein                  |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-60955](http://www.novusbio.com/reviews/submit/NBP2-60955)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

