

Product Datasheet

SOX4 Antibody (CL5634) - BSA Free NBP2-61420

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-61420

Updated 2/27/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-61420



NBP2-61420

SOX4 Antibody (CL5634) - BSA Free


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL5634
Preservative	0.02% Sodium Azide
Isotype	IgG2b
Purity	Protein A purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Host	Mouse
Gene ID	6659
Gene Symbol	SOX4
Species	Human, Mouse, Rat
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions.
Immunogen	This antibody was developed using a synthetic peptide derived from Q06945, with the exact immunogen sequence remaining proprietary.


Product Application Details	
Applications	Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Immunohistochemistry 1:1000 - 1:2500, Immunocytochemistry/ Immunofluorescence 2-10 ug/ml, Immunohistochemistry-Paraffin 1:1000 - 1:2500
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: PFA/Triton X-100.

Images


Immunocytochemistry/Immunofluorescence: SOX4 Antibody (CL5634) [NBP2-61420] - Staining of SH-SY5Y cells using the Anti-SOX4 monoclonal antibody, showing specific staining in the nucleoplasm in green. Microtubule and nuclear probes are visualized in red and blue, respectively (where available). Antibody staining is shown in green.



Immunohistochemistry-Paraffin: SOX4 Antibody (CL5634) [NBP2-61420]
- Staining of rat hippocampus shows strong nuclear immunoreactivity in a subset of cells in the granule cell layer of the dentate gyrus.




Immunohistochemistry-Paraffin: SOX4 Antibody (CL5634) [NBP2-61420]
- Staining of human glioma shows moderate to strong nuclear positivity in tumor cells.




Immunohistochemistry-Paraffin: SOX4 Antibody (CL5634) [NBP2-61420]
- Staining of human liver shows absence of nuclear positivity in hepatocytes as expected (negative control).



Immunohistochemistry-Paraffin: SOX4 Antibody (CL5634) [NBP2-61420]
- Staining of mouse forebrain shows strong positivity in cells in the rostral migratory stream, as well as in a subset of neurons in the olfactory bulb.



Immunohistochemistry-Paraffin: SOX4 Antibody (CL5634) [NBP2-61420]
- Staining of mouse forebrain shows strong immunoreactivity in cells in the rostral migratory stream.



Immunohistochemistry-Paraffin: SOX4 Antibody (CL5634) [NBP2-61420]
- Staining of mouse hippocampus shows strong immunoreactivity in a subset of neurons in the granule cell layer of the dentate gyrus.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-61420

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)
H00006659-Q01-10ug	Recombinant Human SOX4 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-61420

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

