

Product Datasheet

Recombinant TEV Protease His (N-Term) Protein NBP2-61441-300IU

Unit Size: 300IU

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-61441

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-61441



NBP2-61441-300IU**Recombinant TEV Protease His (N-Term) Protein**

Product Information	
Unit Size	300IU
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Preservative	50% glycerol
Purity	>90%, by SDS-PAGE
Buffer	0.2 um filtered solution in 25 mM Tris-HCl, pH 8.0, 75 mM NaCl, 5 mM EDTA, 10 mM GSH
Target Molecular Weight	27 kDa

Product Description	
Description	<p>TEV protease encoded by the tobacco etch virus is a catalytic domain of the Nuclear Inclusion a (Nla) protein. It consists of 241 a.a. amino acids with the molecular weight of 27kDa. It contains N-terminal His-tag.</p> <p>Source: <i>Escherichia coli</i>.</p> <p>One unit is defined as the amount of enzyme needed to cleave 3 ug of fusion protein in 1 hour to 85 % completion at 30C in a buffer containing 50 mM Tris-HCl, pH 8.0, 0.5 mM EDTA, and 1 mM DTT.</p> <p>TEV recognizes the amino acid sequence of the general form E-X-X-Y-X-Q (or S)/X', and cleaves between Q (or S)/X'. In this form X and X' stand for any of the amino acid residues, except that X' cannot be P. The optimal cleavage site is ENLYFQ/G.</p>
Gene ID	1502321
Gene Symbol	TEVgp1
Details of Functionality	TEV recognizes the amino acid sequence of the general form E-X-X-Y-X-Q (or S)/X', and cleaves between Q (or S)/X'. In this form X and X' stand for any of the amino acid residues, except that X' cannot be P. The optimal cleavage site is ENLYFQ/G. As having the absolute specificity and widely using conditions like broad pH range and ionic strength, the TEV protease became more versatile than EK, thrombin and other protease used in biochemical applications, especially recombinant protein production.
Inhibitor Target	Caspase 8

Product Application Details	
Applications	Enzyme Activity
Recommended Dilutions	Enzyme Activity
Application Notes	The optimal temperature for cleavage is 30C; however, the enzyme can be used at temperatures as low as 4C.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-61441-300IU

4469-TP-200	TEV Protease [Unconjugated]
-------------	-----------------------------

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-61441

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

