

Product Datasheet

Histone H3 [p Ser10, Trimethyl Lys9] Antibody (RM162) NBP2-61522

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-61522

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-61522



NBP2-61522

Histone H3 [p Ser10, Trimethyl Lys9] Antibody (RM162)

Product Information	
Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	RM162
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	50% Glycerol/PBS, 1% BSA
Target Molecular Weight	15 kDa


Product Description	
Host	Rabbit
Gene ID	126961
Gene Symbol	H3C14
Species	Human, Vertebrate
Specificity/Sensitivity	This Histone H3 [p Ser10, Trimethyl Lys9] antibody (RM162) reacts to Histone H3 only when modified AT both trimethyl lysine 9 and phospho serine 10.
Immunogen	This Histone H3 [p Ser10, Trimethyl Lys9] antibody (RM162) was raised against a trimethyl-phospho-peptide corresponding to Histone H3 [p Ser10, Trimethyl Lys9]

Product Application Details	
Applications	Western Blot, ELISA, Multiplex Immunoassay
Recommended Dilutions	Western Blot 0.01 ug/ml-1 ug/ml, ELISA 0.01 ug/mL - 0.5 ug/mL, Multiplex Immunoassay 0.1 ug/ml - 1 ug/ml




Images

Western Blot: Histone H3 [p Ser10, Trimethyl Lys9] Antibody (RM162) [NBP2-61522] - Western Blot of acid extracts of HeLa cells (1) and recombinant histone H3.3 (2), using NBP2-61522 at 0.01 ug/ml. Observed molecular weight is ~17 kDa.



Multiplex Immunoassay: Histone H3 [p Ser10, Trimethyl Lys9] Antibody (RM162) [NBP2-61522] - NBP2-61522 specifically reacts to Histone H3 only when modified at both trimethyl-lysine 9 and phospho-serine 10 (K9me3/S10p). No cross reactivity with non-modified Lysine 9/ Serine 10, methylated Lysine 9 (K9me1, k9me2, k9me3) ONLY, or phosphorylated Serine 10 ONLY, in Histone H3.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-61522

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NB21-1253PEP	Histone H3 [Trimethyl Lys36] Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-61522

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

