

Product Datasheet

HIST1H3A Antibody (CL4968) - BSA Free NBP2-61643

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-61643

Updated 2/27/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-61643



NBP2-61643


HIST1H3A Antibody (CL4968) - BSA Free

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL4968
Preservative	0.02% Sodium Azide
Isotype	IgG2b
Purity	Protein A purified
Buffer	PBS (pH 7.2) and 40% Glycerol
Product Description	
Host	Mouse
Gene ID	8350
Gene Symbol	H3C1
Species	Human
Immunogen	This antibody was developed using a synthetic peptide derived from P68431, with the exact immunogen sequence remaining proprietary.
Product Application Details	
Applications	Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Immunohistochemistry 1:200 - 1:500, Immunohistochemistry-Paraffin 1:200 - 1:500
Application Notes	IHC-Paraffin HIER pH6 retrieval is recommended.




Images


Immunohistochemistry: Histone H3T Antibody (CL4968) [NBP2-61643] - Staining of human cerebral cortex shows moderate to strong nuclear immunoreactivity in neuronal and glial cells.




Immunohistochemistry: Histone H3T Antibody (CL4968) [NBP2-61643] - Staining of human prostate shows moderate nuclear positivity in glandular and smooth muscle cells.



Immunohistochemistry: Histone H3T Antibody (CL4968) [NBP2-61643] - Staining of human duodenum shows strong nuclear immunoreactivity in glandular cells as well as in cells in lamina propria.



Immunohistochemistry: Histone H3T Antibody (CL4968) [NBP2-61643] - Staining of human fallopian tube shows strong nuclear positivity in glandular and stromal cells.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-61643

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)
NBP2-06572	HIST2H3C Overexpression Lysate

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-61643

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

