

Product Datasheet

LAIR1 Antibody (NKTA255) - BSA Free NBP2-62214

Unit Size: 0.1 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Reviews: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-62214

Updated 2/23/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-62214



NBP2-62214


LAIR1 Antibody (NKTA255) - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	NKTA255
Preservative	0.09% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH 7.4)
Target Molecular Weight	40 kDa
Product Description	
Host	Mouse
Gene ID	3903
Gene Symbol	LAIR1
Species	Human
Specificity/Sensitivity	The mouse monoclonal antibody NKTA255 recognizes CD305 / LAIR1, a 40 kDa type I transmembrane glycoprotein expressed on NK, T, and B cells, monocytes, dendritic cells, eosinophils, basophils, mast cells, CD34+ hematopoietic progenitor cells and thymocytes.
Immunogen	Activated NK cells and CD3- thymocytes
Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:100-1:2000, Flow Cytometry 1-4 ug/ml, Immunoprecipitation 1:10-1:500




Images


Flow (Cell Surface): LAIR1 Antibody (NKTA255) [NBP2-62214] - Surface staining of CD305 / LAIR-1 in human peripheral blood with anti-CD305 (NKTA255) PerCP-CyTM5.5.



Flow (Cell Surface): LAIR1 Antibody (NKTA255) [NBP2-62214] - Surface staining of CD305 / LAIR-1 in human peripheral blood (lymphocyte gate) with anti-CD305 (NKTA255) PE.



Surface staining of human lymphocytes stained using anti-human CD305 (NKTA255) purified antibody (concentration in sample 2 ug/ml, GAM APC) and anti-human CD19 (LT19) PE antibody (20 µl reagent / 100 µl of peripheral whole blood).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-62214

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP1-84590PEP	LAIR1 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-62214

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

