

Product Datasheet

Human alpha-Synuclein ELISA Kit (Chemiluminescence) NBP2-62771

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-62771

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-62771




NBP2-62771**Human alpha-Synuclein ELISA Kit (Chemiluminescence)**


Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.
Product Description	
Gene ID	6622
Gene Symbol	SNCA
Species	Human
Specificity/Sensitivity	This Human alpha-Synuclein ELISA Kit (Chemiluminescence) is specific for Human alpha-Synuclein in samples. No significant cross-reactivity or interference between Human alpha-Synuclein and analogues was observed.
Kit Components	Micro CLIA Plate, Reference Standard, Concentrated Biotinylated Detection Ab (100x), Concentrated HRP Conjugate (100x), Sample Diluent, Biotinylated Detection Ab Diluent, HRP Conjugate Diluent, Concentrated Wash Buffer (25x), Substrate Reagent A, Substrate Reagent B, Plate Sealer, Product Manual
Standard Curve Range	12.50-800 pg/mL
Sensitivity	7.50 pg/mL
Inter Assay	CV% < 10.2%
Intra Assay	CV% < 10.2%
Assay Type	Sandwich-CLIA
Spiking Recovery	86-108%
Suitable Sample Type	Serum, plasma and other biological fluids
Sample Volume	100 uL
Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA

Images

ELISA: Human alpha-Synuclein ELISA Kit (Chemiluminescence) [NBP2-62771] - Samples were spiked with high concentrations of Human alpha-Synuclein and diluted with Reference Standard & Sample Diluent to produce samples with values within the range of the assay.



ELISA: Human alpha-Synuclein ELISA Kit (Chemiluminescence) [NBP2-62771] - Standard Curve Reference





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-62771

NBP2-54789-100ug	Recombinant Human alpha-Synuclein Active, Pre-formed Fibrils, (Type 1) Protein
BC100-494	PINK1 Antibody - BSA Free
NBP2-15365	alpha-Synuclein Antibody
NB300-109	Tyrosine Hydroxylase Antibody

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-62771

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

