

Product Datasheet

Cytokeratin 18 Antibody (SZ80-07) NBP2-67370

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-67370

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-67370



NBP2-67370**Cytokeratin 18 Antibody (SZ80-07)**


| Product Information | |
|---------------------|--|
| Unit Size | 100 ul |
| Concentration | 1 mg/ml |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | SZ80-07 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG |
| Purity | Protein A purified |
| Buffer | TBS (pH7.4), 0.05% BSA, 40% Glycerol |

| Product Description | |
|---------------------|---|
| Host | Rabbit |
| Gene ID | 3875 |
| Gene Symbol | KRT18 |
| Species | Human, Mouse, Rat |
| Immunogen | Recombinant protein within mouse Cytokeratin 18 aa 200-423. (SwissProt: P05783 Human; SwissProt: P05784 Mouse; SwissProt: Q5BJY9 Rat) |


| Product Application Details | |
|-----------------------------|---|
| Applications | Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin |
| Recommended Dilutions | Western Blot 1:500-1:2000, Flow Cytometry 1:50-1:100, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1:100-1:500, Immunohistochemistry-Paraffin 1:50-1:200 |

Images


Western Blot: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Western blot analysis of Cytokeratin 18 on different lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody (1/500) was used in 5% BSA at room temperature for 2 hours. Goat




Immunocytochemistry/Immunofluorescence: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Staining Cytokeratin 18 in HepG2 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Immunohistochemistry-Paraffin: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Analysis of paraffin-embedded human liver tissue with Rabbit anti-Cytokeratin 18 antibody washed with ddH2O and PBS, and then probed with the primary antibody (1/200) for 30 minutes at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.




Flow Cytometry: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Analysis of MCF-7 cells with Cytokeratin 18 antibody at 1/50 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti rabbit IgG was used as the secondary antibody




Immunocytochemistry/Immunofluorescence: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Staining Cytokeratin 18 in A431 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Immunocytochemistry/Immunofluorescence: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Staining Cytokeratin 18 in Hela cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



Immunohistochemistry-Paraffin: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Analysis of paraffin-embedded human breast carcinoma tissue using anti-Cytokeratin 18 antibody. Counter stained with hematoxylin.




Immunohistochemistry-Paraffin: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Analysis of paraffin-embedded human colon cancer tissue using anti-Cytokeratin 18 antibody. Counter stained with hematoxylin.




Immunohistochemistry-Paraffin: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Analysis of paraffin-embedded human gastric carcinoma tissue using anti-Cytokeratin 18 antibody. Counter stained with hematoxylin.




Immunohistochemistry-Paraffin: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Analysis of paraffin-embedded human kidney tissue using anti-Cytokeratin 18 antibody. Counter stained with hematoxylin.




Immunohistochemistry-Paraffin: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Analysis of paraffin-embedded human liver tissue using anti-Cytokeratin 18 antibody. Counter stained with hematoxylin.




Immunohistochemistry-Paraffin: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Analysis of paraffin-embedded human colon tissue with Rabbit anti-Cytokeratin 18 antibody washed with ddH₂O and PBS, and then probed with the primary antibody (1/200) for 30 minutes at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.




Immunohistochemistry-Paraffin: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Analysis of paraffin-embedded human stomach carcinoma tissue with Rabbit anti-Cytokeratin 18 antibody washed with ddH₂O and PBS, and then probed with the primary antibody (1/200) for 30 minutes at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.



Immunohistochemistry-Paraffin: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Analysis of paraffin-embedded human stomach carcinoma tissue with Rabbit anti-Cytokeratin 18 antibody washed with ddH₂O and PBS, and then probed with the primary antibody (1/200) for 30 minutes at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.



Immunohistochemistry-Paraffin: Cytokeratin 18 Antibody (SZ80-07) [NBP2-67370] - Analysis of paraffin-embedded human colon carcinoma with Rabbit anti-Cytokeratin 18 antibody washed with ddH₂O and PBS, and then probed with the primary antibody (1/200) for 30 minutes at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.



Publications

Lodes NT, Seidensticker K, Perniss A et al. Investigation on ciliary functionality of different airway epithelial cell lines in 3D cell culture Tissue Eng Part A 2019-11-07 [PMID: 31696788] (IHC-P, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-67370

| | |
|--------------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| H00003875-P02-10ug | Recombinant Human Cytokeratin 18 GST (N-Term) Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-67370

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

