

Product Datasheet

B7-2/CD86 Antibody (SJ20-00) NBP2-67417

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-67417

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-67417



NBP2-67417


B7-2/CD86 Antibody (SJ20-00)

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	SJ20-00
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	TBS (pH7.4), 0.05% BSA, 40% Glycerol
Product Description	
Host	Rabbit
Gene ID	942
Gene Symbol	CD86
Species	Human, Mouse, Rat
Immunogen	Synthetic peptide within Human B7-2/CD86 aa 1-50 / 329. (SwissProt: P42081 Human; Entrez Gene:56822 Rat)
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:2000, Immunocytochemistry/ Immunofluorescence 1:500, Immunoprecipitation Use at an assay dependent concentration




Images


Immunocytochemistry/Immunofluorescence: B7-2/CD86 Antibody (SJ20-00) [NBP2-67417] - Analysis of Raji cells labeling CD86 with Rabbit anti-CD86 antibody at 1/500 dilution. Cells were fixed in 100% precooled methanol for 5 minutes at room temperature, then blocked with 1% BSA in 10% negative goat serum for 1 hour at room temperature. Cells were then incubated with Rabbit anti-CD86 antibody at 1/500 dilution in 1% BSA in PBST overnight at 4 °C. Goat Anti-Rabbit IgG H&L (iFluor™ 488, was used as the secondary antibody at 1/1,000 dilution. PBS instead of the primary antibody was used as the secondary antibody only control. Nuclear DNA was labelled in blue with DAPI. Beta tubulin (red) was stained at 1/100 dilution overnight at +4 °C. Goat Anti-Mouse IgG H&L (iFluor™ 594, was used as the secondary antibody at 1/1,000 dilution.



Western Blot: B7-2/CD86 Antibody (SJ20-00) [NBP2-67417] - Analysis of CD86 on different lysates with Rabbit anti-CD86 antibody at 1/1,000 dilution. Lane 1: Raji cell lysate (10 ug/Lane) Lane 2: Mouse lung tissue lysate (20 ug/Lane) Lane 3: Rat lung tissue lysate (20 ug/Lane)
 Predicted band size: 38 kDa Observed band size: 70 kDa Exposure time: 1 minute; 4-20% SDS-PAGE gel. Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody at 1/1,000 dilution was used in 5% NFDM/TBST at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody at 1:100,000 dilution was used for 1 hour at room temperature.



Western Blot: B7-2/CD86 Antibody (SJ20-00) [NBP2-67417] - Analysis of CD86 on different lysates with Rabbit anti-CD86 antibody at 1/2,000 dilution. Lane 1: Daudi cell lysate Lane 2: Raji cell lysate Lane 3: Ramos cell lysate Lane 4: Jurkat cell lysate Lane 5: K-562 cell lysate
 Lysates/proteins at 15 ug/Lane. Predicted band size: 38 kDa Observed band size: 70 kDa Exposure time: 2 minutes 37 seconds; 4-20% SDS-PAGE gel. Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody at 1/2,000 dilution was used in 5% NFDM/TBST at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody at 1:100,000 dilution was used for 1 hour at room temperature.



Publications

Kadowaki M, Yoshida S, Itoyama T et al. Involvement of M1/M2 Macrophage Polarization in Reparative Dentin Formation Life (Basel, Switzerland) 2022-11-07 [PMID: 36362965] (IHC-P, Rat)

Details:
 Dilution used for IHC-P 1:50

Zavatti M, Beretti F, Casciaro F et al. Comparison of the therapeutic effect of amniotic fluid stem cells and their exosomes on monoiodoacetate-induced animal model of osteoarthritis Biofactors 2019-10-18 [PMID: 31625201] (IF/IHC)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-67417

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-77486-50ug	Recombinant Human B7-2/CD86 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-67417

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

