

# Product Datasheet

## Goat anti-Mouse IgG3 Secondary Antibody [PerCP] NBP2-68520PCP

Unit Size: 0.5 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-68520PCP](http://www.novusbio.com/NBP2-68520PCP)

Updated 9/6/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-68520PCP](http://www.novusbio.com/reviews/destination/NBP2-68520PCP)



**NBP2-68520PCP**

Goat anti-Mouse IgG3 Secondary Antibody [PerCP]

| Product Information         |                                                                                                                                                                                                                                                         |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.5 ml                                                                                                                                                                                                                                                  |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                             |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                              |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                      |
| Isotype                     | IgG                                                                                                                                                                                                                                                     |
| Conjugate                   | PerCP                                                                                                                                                                                                                                                   |
| Purity                      | Affinity purified                                                                                                                                                                                                                                       |
| Buffer                      | PBS                                                                                                                                                                                                                                                     |
| Product Description         |                                                                                                                                                                                                                                                         |
| Host                        | Goat                                                                                                                                                                                                                                                    |
| Gene ID                     | 3502                                                                                                                                                                                                                                                    |
| Gene Symbol                 | IGHG3                                                                                                                                                                                                                                                   |
| Species                     | Mouse                                                                                                                                                                                                                                                   |
| Specificity/Sensitivity     | Goat anti Mouse IgG3 antibody recognizes Mouse IgG3 and has been pre-absorbed against mouse IgM, IgG1, IgG2a, IgG2b and IgA, pooled human sera and purified human paraproteins. The antibody shows minimal cross-reactivity with human immunoglobulins. |
| Immunogen                   | Mouse IgG3 paraproteins                                                                                                                                                                                                                                 |
| Product Application Details |                                                                                                                                                                                                                                                         |
| Applications                | ELISA                                                                                                                                                                                                                                                   |
| Recommended Dilutions       | ELISA                                                                                                                                                                                                                                                   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                  |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP2-68520PCP**

---

|                    |                                                    |
|--------------------|----------------------------------------------------|
| H00003502-P01-10ug | Recombinant Human IgG3 GST (N-Term) Protein        |
| 210-TA-005         | TNF-alpha [Unconjugated]                           |
| MAB007             | Mouse IgG3 Isotype Control (133316) [Unconjugated] |
| M6000B-1           | IL-6 [HRP]                                         |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Secondary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-68520PCP](http://www.novusbio.com/reviews/submit/NBP2-68520PCP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

