

# Product Datasheet

## Dolichol Kinase Antibody - BSA Free NBP2-68873

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-68873](http://www.novusbio.com/NBP2-68873)

Updated 2/23/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-68873](http://www.novusbio.com/reviews/destination/NBP2-68873)



**NBP2-68873**

Dolichol Kinase Antibody - BSA Free

| Product Information |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 100 ul                                                                                                           |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                           |
| Clonality           | Polyclonal                                                                                                       |
| Preservative        | 0.02% Sodium Azide                                                                                               |
| Isotype             | IgG                                                                                                              |
| Purity              | Protein A purified                                                                                               |
| Buffer              | PBS (pH 7.2) and 40% Glycerol                                                                                    |


| Product Description |                                                                                                                 |
|---------------------|-----------------------------------------------------------------------------------------------------------------|
| Host                | Rabbit                                                                                                          |
| Gene ID             | 22845                                                                                                           |
| Gene Symbol         | DOLK                                                                                                            |
| Species             | Human                                                                                                           |
| Immunogen           | This antibody was developed against a recombinant protein corresponding to amino acids: YQNAKRSSSESKKHQAPTIARKY |

| Product Application Details |                                                                                   |
|-----------------------------|-----------------------------------------------------------------------------------|
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence                             |
| Recommended Dilutions       | Western Blot 0.04-0.4 ug/ml, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml |
| Application Notes           | ICC/IF, Fixation/Permeabilization: PFA/Triton X-100                               |




## Images

Western Blot: Dolichol Kinase Antibody [NBP2-68873] - Analysis in human cell line RT-4 and human cell line U-251 MG.



Immunocytochemistry/Immunofluorescence: Dolichol Kinase Antibody [NBP2-68873] - Staining of human cell line HaCaT shows localization to nucleus, nucleoli & cytosol.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-68873**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| NBP2-68873PEP | Dolichol Kinase Recombinant Protein Antigen         |
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-68873](http://www.novusbio.com/reviews/submit/NBP2-68873)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

