

Product Datasheet

Goat anti-Mouse IgD Fc Secondary Antibody NBP2-69334

Unit Size: 10 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-69334

Updated 5/14/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-69334



NBP2-69334**Goat anti-Mouse IgD Fc Secondary Antibody**

Product Information	
Unit Size	10 mg
Concentration	LYOPH mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Reconstitution Instructions	Reconstitute by adding 1 ml sterile distilled water, spun down to remove insoluble particles, divided into small aliquots, frozen and stored at or below -20C.
Purity	Protein G purified
Buffer	Lyophilized from PBS (pH 7.2)

Product Description	
Description	The lyophilized 7S (IgG) fraction is shipped at ambient temperature and may be stored at 4C; prolonged storage at or below -20C. Prior to use, an aliquot is thawed slowly at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at 4C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the product.
Host	Goat
Species	Mouse
Reactivity Notes	Inter-species cross-reactivity is a normal feature of antibodies to immunoglobulins, since Ig of different species frequently share antigenic determinants. Cross-reactivity of this immunoconjugate has not been tested in detail.
Specificity/Sensitivity	The reactivity of the antiserum is directed to the Fc subunit of the IgD molecule which expresses strict isotypic (class) specificity. In immunoelectrophoresis and radial immunodiffusion, using various antiserum concentrations against serum of mice belonging to different allotypic groups, a single precipitin line is obtained which show a reaction of identity with the precipitin lines obtained with the purified IgD of BALB/C and C57BL origin used as immunogens. It does not react with IgG, including all subclasses, IgM and IgG/Fab or any non- Ig protein in mouse serum, as tested by immunoelectrophoresis and double radial immunodiffusion. The antiserum also reacts with membrane-bound IgD in peripheral blood cells of different mouse strains as tested by immunofluorescence microscopy.
Immunogen	Purified polyclonal and monoclonal IgD isolated from BALB/C and C57BL mouse serum. Freund's complete adjuvant is used in the first step of the immunization procedure.

Product Application Details	
Applications	Western Blot, Dot Blot, ELISA, Immunoblotting, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot, ELISA 1:100-1:2000, Immunocytochemistry/ Immunofluorescence 1:10 - 1:500, Immunoblotting 1:100 - 1:2000, Dot Blot 1:100-1:2000
Application Notes	Indirect immunofluorescence.

Publications

Yang Zhang, Laura Garcia-Ibanez, Carolin Ulbricht, Laurence S. C. Lok, Jeremy A. Pike, Jennifer Mueller-Winkler, Thomas W. Dennison, John R. Ferdinand, Cameron J. M. Burnett, Juan C. Yam-Puc, Lingling Zhang, Raul Maqueda Alfaro, Yousuke Takahama, Izumi Ohigashi, Geoffrey Brown, Tomohiro Kurosaki, Victor L. J. Tybulewicz, Antal Rot, Anja E. Hauser, Menna R. Clatworthy, Kai-Michael Toellner Recycling of memory B cells between germinal center and lymph node subcapsular sinus supports affinity maturation to antigenic drift Nature Communications 2022-05-05 [PMID: 35513371]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Secondary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-69334

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

