

# Product Datasheet

## Catalase Antibody (OT11B8) [DyLight 405] NBP2-70342V

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-70342V](http://www.novusbio.com/NBP2-70342V)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-70342V](http://www.novusbio.com/reviews/destination/NBP2-70342V)



**NBP2-70342V**

Catalase Antibody (OTI1B8) [DyLight 405]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                 |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                          |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                     |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                        |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                      |
| Clone                       | OTI1B8                                                                                                                                                                                                                                                                                                                                          |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                              |
| Isotype                     | IgG1                                                                                                                                                                                                                                                                                                                                            |
| Conjugate                   | DyLight 405                                                                                                                                                                                                                                                                                                                                     |
| Purity                      | Immunogen affinity purified                                                                                                                                                                                                                                                                                                                     |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                              |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                 |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                           |
| Gene ID                     | 847                                                                                                                                                                                                                                                                                                                                             |
| Gene Symbol                 | CAT                                                                                                                                                                                                                                                                                                                                             |
| Species                     | Human, Mouse, Rat, Monkey, Primate                                                                                                                                                                                                                                                                                                              |
| Reactivity Notes            | Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions. |
| Immunogen                   | Full length recombinant protein of human Catalase (NP_111743) produced in HEK293T cells.                                                                                                                                                                                                                                                        |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                                               |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                 |
| Applications                | Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                            |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                            |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                          |



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-70342V**

---

|                   |                                                |
|-------------------|------------------------------------------------|
| NBP1-97005V-0.5ml | Mouse IgG1 Isotype Control (MG1) [DyLight 405] |
| NBP1-98907-100ug  | Recombinant Human Catalase His Protein         |
| 210-TA-005        | TNF-alpha [Unconjugated]                       |
| NBP3-24476        | Catalase Activity Assay Kit (Colorimetric)     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-70342V](http://www.novusbio.com/reviews/submit/NBP2-70342V)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

