

# Product Datasheet

## CRY2 Antibody (OT1H5) - Azide and BSA Free NBP2-70449

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-70449](http://www.novusbio.com/NBP2-70449)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-70449](http://www.novusbio.com/reviews/destination/NBP2-70449)



**NBP2-70449**


CRY2 Antibody (OTI1H5) - Azide and BSA Free

| Product Information         |                                                                                                                                                                                                                                                                                                                                                 |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ug                                                                                                                                                                                                                                                                                                                                          |
| Concentration               | LYOPH mg/ml                                                                                                                                                                                                                                                                                                                                     |
| Storage                     | Store at -20C. Avoid freeze-thaw cycles.                                                                                                                                                                                                                                                                                                        |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                      |
| Clone                       | OTI1H5                                                                                                                                                                                                                                                                                                                                          |
| Preservative                | No Preservative                                                                                                                                                                                                                                                                                                                                 |
| Reconstitution Instructions | we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.                                                                                                                |
| Isotype                     | IgG2a                                                                                                                                                                                                                                                                                                                                           |
| Purity                      | Immunogen affinity purified                                                                                                                                                                                                                                                                                                                     |
| Buffer                      | Lyophilized from PBS (pH 7.3) with 8% Trehalose                                                                                                                                                                                                                                                                                                 |
| Target Molecular Weight     | 66.8 kDa                                                                                                                                                                                                                                                                                                                                        |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                 |
| Description                 | Novus Biologicals Mouse CRY2 Antibody (OTI1H5) - Azide and BSA Free (NBP2-00749) is a monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                                                                                    |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                           |
| Gene ID                     | 1408                                                                                                                                                                                                                                                                                                                                            |
| Gene Symbol                 | CRY2                                                                                                                                                                                                                                                                                                                                            |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                               |
| Reactivity Notes            | Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions. |
| Immunogen                   | Full length human recombinant protein of human CRY2 (NP_0066940) produced in HEK293T cell.                                                                                                                                                                                                                                                      |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                 |
| Applications                | Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, CyTOF-ready                                                                                                                                                                                                         |
| Recommended Dilutions       | Western Blot 1:2000, Flow Cytometry 1:100, Immunohistochemistry 1:150, Immunocytochemistry/ Immunofluorescence 1:100, Immunohistochemistry-Paraffin, CyTOF-ready                                                                                                                                                                                |




## Images


Western Blot: CRY2 Antibody (OTI1H5) - Azide and BSA Free [NBP2-70449] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CRY2 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CRY2.




Immunocytochemistry/Immunofluorescence: CRY2 Antibody (OTI1H5) - Azide and BSA Free [NBP2-70449] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY CRY2.




Immunohistochemistry: CRY2 Antibody (OTI1H5) - Azide and BSA Free [NBP2-70449] - Staining of paraffin-embedded Human colon tissue using anti-CRY2 mouse monoclonal antibody.



Flow Cytometry: CRY2 Antibody (OTI1H5) - Azide and BSA Free [NBP2-70449] - Analysis of Jurkat cells, using anti-CRY2 antibody, (Red), compared to a nonspecific negative control antibody (Blue).




Immunohistochemistry: CRY2 Antibody (OTI1H5) - Azide and BSA Free [NBP2-70449] - Staining of paraffin-embedded Adenocarcinoma of Human breast tissue using anti-CRY2 mouse monoclonal antibody.




Immunohistochemistry: CRY2 Antibody (OTI1H5) - Azide and BSA Free [NBP2-70449] - Staining of paraffin-embedded Adenocarcinoma of Human colon tissue using anti-CRY2 mouse monoclonal antibody.



Immunohistochemistry: CRY2 Antibody (OTI1H5) - Azide and BSA Free [NBP2-70449] - Staining of paraffin-embedded Carcinoma of Human bladder tissue using anti-CRY2 mouse monoclonal antibody.




Immunohistochemistry: CRY2 Antibody (OTI1H5) - Azide and BSA Free [NBP2-70449] - Staining of paraffin-embedded Carcinoma of Human pancreas tissue using anti-CRY2 mouse monoclonal antibody.






Immunohistochemistry: CRY2 Antibody (OTI1H5) - Azide and BSA Free [NBP2-70449] - Staining of paraffin-embedded Carcinoma of Human thyroid tissue using anti-CRY2 mouse monoclonal antibody.



Immunohistochemistry: CRY2 Antibody (OTI1H5) - Azide and BSA Free [NBP2-70449] - Staining of paraffin-embedded Human lymphoma tissue using anti-CRY2 mouse monoclonal antibody.



Flow Cytometry: CRY2 Antibody (OTI1H5) - Azide and BSA Free [NBP2-70449] - Analysis of Hela cells, using anti-CRY2 antibody, (Red), compared to a nonspecific negative control antibody (Blue).





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-70449**

|                    |                                                         |
|--------------------|---------------------------------------------------------|
| HAF007             | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B            | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-96778         | Mouse IgG2a Isotype Control (M2A)                       |
| H00001408-Q01-10ug | Recombinant Human CRY2 GST (N-Term) Protein             |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-70449](http://www.novusbio.com/reviews/submit/NBP2-70449)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

