

# Product Datasheet

## Cytochrome P450 2B6 Antibody (OTI3D5) [DyLight 680] NBP2-70526FR

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-70526FR](http://www.novusbio.com/NBP2-70526FR)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-70526FR](http://www.novusbio.com/reviews/destination/NBP2-70526FR)



**NBP2-70526FR**

Cytochrome P450 2B6 Antibody (OTI3D5) [DyLight 680]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                 |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                          |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                     |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                        |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                      |
| Clone                       | OTI3D5                                                                                                                                                                                                                                                                                                                                          |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                              |
| Isotype                     | IgG2a                                                                                                                                                                                                                                                                                                                                           |
| Conjugate                   | DyLight 680                                                                                                                                                                                                                                                                                                                                     |
| Purity                      | Immunogen affinity purified                                                                                                                                                                                                                                                                                                                     |
| Buffer                      | PBS                                                                                                                                                                                                                                                                                                                                             |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                 |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                           |
| Gene ID                     | 1555                                                                                                                                                                                                                                                                                                                                            |
| Gene Symbol                 | CYP2B6                                                                                                                                                                                                                                                                                                                                          |
| Species                     | Human, Mouse, Rat, Canine, Monkey                                                                                                                                                                                                                                                                                                               |
| Reactivity Notes            | Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions. |
| Immunogen                   | Full length recombinant protein of human Cytochrome P450 2B6 (NP_000758) produced in HEK293T cells.                                                                                                                                                                                                                                             |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                                               |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                 |
| Applications                | Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                            |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                            |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                          |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-70526FR**

---

|                    |                                                            |
|--------------------|------------------------------------------------------------|
| H00001555-P01-10ug | Recombinant Human Cytochrome P450 2B6 GST (N-Term) Protein |
| NBP2-04177         | Cytochrome P450 2B6 Overexpression Lysate                  |
| NB120-19347        | Calcium-sensing R/CaSR Antibody (5C10, ADD) - BSA Free     |
| NBP1-85367         | Cytochrome P450 2E1 Antibody                               |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-70526FR](http://www.novusbio.com/reviews/submit/NBP2-70526FR)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

