

Product Datasheet

Cytokeratin 18 Antibody (OT1E1) [mFluor Violet 500 SE] NBP2-71086MFV500

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-71086MFV500

Updated 11/9/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-71086MFV500



NBP2-71086MFV500

Cytokeratin 18 Antibody (OT11E1) [mFluor Violet 500 SE]

Product Information

| | |
|----------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | OT11E1 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 |
| Conjugate | mFluor Violet 500 SE |
| Purity | Immunogen affinity purified |
| Buffer | 50mM Sodium Borate |

Product Description

| | |
|-------------------------|---|
| Host | Mouse |
| Gene ID | 3875 |
| Gene Symbol | KRT18 |
| Species | Human, Mouse, Rat, Monkey |
| Reactivity Notes | Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions. |
| Immunogen | Human recombinant protein fragment corresponding to amino acids 69-372 of human CK18 (NP_000215) produced in E.coli. |

Product Application Details

| | |
|------------------------------|--|
| Applications | Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry |
| Recommended Dilutions | Western Blot, Immunohistochemistry, Immunocytochemistry/Immunofluorescence |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |

Images

Cytokeratin 18 Antibody (OT11E1) [mFluor Violet 500 SE] - Vial of mFluor Violet 500 conjugated antibody. mFluor Violet 500 is optimally excited at 410 nm by the Violet laser (405 nm) and has an emission maximum of 501 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-71086MFV500

| | |
|--------------------|---|
| H00003875-P01-10ug | Recombinant Human Cytokeratin 18 GST (N-Term) Protein |
| D6050 | IL-6 [HRP] |
| H00003875-P02-10ug | Recombinant Human Cytokeratin 18 GST (N-Term) Protein |
| AF835 | Caspase-3 Antibody [Unconjugated] - Active |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-71086MFV500

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

