

# Product Datasheet

## Galectin-3 Antibody (OT11C7) [PE] NBP2-71124PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-71124PE](http://www.novusbio.com/NBP2-71124PE)

Updated 1/21/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-71124PE](http://www.novusbio.com/reviews/destination/NBP2-71124PE)



**NBP2-71124PE**

Galectin-3 Antibody (OT11C7) [PE]

| Product Information |   |
|---------------------|---|
| Unit Size           | 0.1 ml  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.  |
| Clonality           | Monoclonal  |
| Clone               | OT11C7  |
| Preservative        | 0.05% Sodium Azide  |
| Isotype             | IgG   |
| Conjugate           | PE  |
| Purity              | Immunogen affinity purified   |
| Buffer              | PBS   |

| Product Description |   |
|---------------------|---|
| Host                | Mouse   |
| Gene ID             | 3958  |
| Gene Symbol         | LGALS3  |
| Species             | Human, Mouse  |
| Reactivity Notes    | Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions. |
| Immunogen           | Full length human recombinant protein of human LGALS3 (NP_002297) produced in E.coli.   |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | Western Blot, Immunohistochemistry                                     |
| Recommended Dilutions       | Western Blot, Immunohistochemistry                                     |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-71124PE**

---

|                 |                                      |
|-----------------|--------------------------------------|
| NBP2-34881-10ug | Recombinant Human Galectin-3 Protein |
| 210-TA-005      | TNF-alpha [Unconjugated]             |
| DGAL30          | Galectin-3 [HRP]                     |
| 6507-IL-010/CF  | IL-4 [Unconjugated]                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-71124PE](http://www.novusbio.com/reviews/submit/NBP2-71124PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

