

# Product Datasheet

## GEMIN2 Antibody (OTI4C11) [FITC] NBP2-71374F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-71374F](http://www.novusbio.com/NBP2-71374F)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-71374F](http://www.novusbio.com/reviews/destination/NBP2-71374F)



**NBP2-71374F**

GEMIN2 Antibody (OTI4C11) [FITC]

| Product Information         |   |
|-----------------------------|---|
| Unit Size                   | 0.1 ml  |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.   |
| Storage                     | Store at 4C in the dark.  |
| Clonality                   | Monoclonal  |
| Clone                       | OTI4C11   |
| Preservative                | 0.05% Sodium Azide  |
| Isotype                     | IgG1  |
| Conjugate                   | FITC  |
| Purity                      | Immunogen affinity purified   |
| Buffer                      | PBS   |
| Product Description         |   |
| Host                        | Mouse   |
| Gene ID                     | 8487  |
| Gene Symbol                 | GEMIN2  |
| Species                     | Human, Mouse, Rat   |
| Reactivity Notes            | Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions. |
| Immunogen                   | Full length recombinant protein of human GEMIN2 (NP_003607) produced in E.coli.   |
| Product Application Details |   |
| Applications                | Western Blot, Immunohistochemistry  |
| Recommended Dilutions       | Western Blot, Immunohistochemistry  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.  |



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-71374F**

---

|               |   |
|---------------|---|
| NBP1-96849    | Mouse IgG1 Isotype Control (MG1) [FITC] |
| NBP2-38553PEP | GEMIN2 Recombinant Protein Antigen      |
| NBL1-15969    | GEMIN2 Overexpression Lysate            |
| AF748         | E-Cadherin Antibody [Unconjugated]      |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-71374F](http://www.novusbio.com/reviews/submit/NBP2-71374F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

