

Product Datasheet

PDLIM2 Antibody (OTI11A12) - Azide and BSA Free NBP2-73311

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-73311

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-73311



NBP2-73311

PDLIM2 Antibody (OTI11A12) - Azide and BSA Free


Product Information	
Unit Size	100 ug
Concentration	LYOPH mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI11A12
Preservative	No Preservative
Reconstitution Instructions	we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Purity	Immunogen affinity purified
Buffer	Lyophilized from PBS (pH 7.3) with 8% Trehalose
Target Molecular Weight	37.3 kDa

Product Description	
Description	Novus Biologicals Mouse PDLIM2 Antibody (OTI11A12) - Azide and BSA Free (NBP2-00619) is a monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	64236
Gene Symbol	PDLIM2
Species	Human
Immunogen	Full length human recombinant protein of human PDLIM2 (NP_067643) produced in HEK293T cell.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, CyTOF-ready
Recommended Dilutions	Western Blot 1:2000, Flow Cytometry 1:100, Immunohistochemistry 1:100, Immunocytochemistry/ Immunofluorescence 1:100, Immunohistochemistry-Paraffin, CyTOF-ready

Images


Western Blot: PDLIM2 Antibody (OTI11A12) - Azide and BSA Free [NBP2-73311] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PDLIM2 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PDLIM2.




Immunocytochemistry/Immunofluorescence: PDLIM2 Antibody (OTI11A12) - Azide and BSA Free [NBP2-73311] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY PDLIM2.




Immunohistochemistry: PDLIM2 Antibody (OTI11A12) - Azide and BSA Free [NBP2-73311] - Staining of paraffin-embedded Human lymph node tissue using anti-PDLIM2 mouse monoclonal antibody.




Flow Cytometry: PDLIM2 Antibody (OTI11A12) - Azide and BSA Free [NBP2-73311] - HEK293T cells transfected with either overexpression plasmid (Red) or empty vector control plasmid (Blue) were immunostaining by anti-PDLIM2 antibody, and then analyzed by flow cytometry.




Immunohistochemistry: PDLIM2 Antibody (OTI11A12) - Azide and BSA Free [NBP2-73311] - Staining of paraffin-embedded Human Kidney tissue using anti-PDLIM2 mouse monoclonal antibody.




Immunohistochemistry: PDLIM2 Antibody (OTI11A12) - Azide and BSA Free [NBP2-73311] - Staining of paraffin-embedded Human liver tissue using anti-PDLIM2 mouse monoclonal antibody.




Immunohistochemistry: PDLIM2 Antibody (OTI11A12) - Azide and BSA Free [NBP2-73311] - Staining of paraffin-embedded Human lymphoma tissue using anti-PDLIM2 mouse monoclonal antibody.



Flow Cytometry: PDLIM2 Antibody (OTI11A12) - Azide and BSA Free [NBP2-73311] - Analysis of Hela cells, using anti-PDLIM2 antibody, (Red), compared to a nonspecific negative control antibody (Blue).



Flow Cytometry: PDLIM2 Antibody (OTI11A12) - Azide and BSA Free [NBP2-73311] - Analysis of Jurkat cells, using anti-PDLIM2 antibody, (Red), compared to a nonspecific negative control antibody (Blue).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-73311

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-73311

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

