

# Product Datasheet

## RPA2 Antibody (OTI9A1) - Azide and BSA Free NBP2-73935

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-73935](http://www.novusbio.com/NBP2-73935)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-73935](http://www.novusbio.com/reviews/destination/NBP2-73935)



**NBP2-73935**


RPA2 Antibody (OTI9A1) - Azide and BSA Free

<b>Product Information</b>	
<b>Unit Size</b>	100 ug
<b>Concentration</b>	LYOPH mg/ml
<b>Storage</b>	Store at -20C. Avoid freeze-thaw cycles.
<b>Clonality</b>	Monoclonal
<b>Clone</b>	OTI9A1
<b>Preservative</b>	No Preservative
<b>Reconstitution Instructions</b>	we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
<b>Isotype</b>	IgG1
<b>Purity</b>	Immunogen affinity purified
<b>Buffer</b>	Lyophilized from PBS (pH 7.3) with 8% Trehalose
<b>Target Molecular Weight</b>	29.1 kDa
<b>Product Description</b>	
<b>Description</b>	Novus Biologicals Mouse RPA2 Antibody (OTI9A1) - Azide and BSA Free (NBP2-02182) is a monoclonal antibody validated for use in IHC, WB, Flow, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
<b>Host</b>	Mouse
<b>Gene ID</b>	6118
<b>Gene Symbol</b>	RPA2
<b>Species</b>	Human, Mouse, Rat, Monkey, Primate
<b>Reactivity Notes</b>	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions.
<b>Immunogen</b>	Full length human recombinant protein of human RPA2 (NP_002937) produced in HEK293T cell.
<b>Product Application Details</b>	
<b>Applications</b>	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation, CyTOF-ready
<b>Recommended Dilutions</b>	Western Blot 1:1000-2000, Flow Cytometry 1:100, Immunohistochemistry 1:50, Immunocytochemistry/ Immunofluorescence 1:100, Immunoprecipitation 2ug/500ul, Immunohistochemistry-Paraffin, CyTOF-ready




## Images


Western Blot: RPA2 Antibody (OTI9A1) - Azide and BSA Free [NBP2-73935] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RPA2 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RPA2.




Immunocytochemistry/Immunofluorescence: RPA2 Antibody (OTI9A1) - Azide and BSA Free [NBP2-73935] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY RPA2.




Immunohistochemistry: RPA2 Antibody (OTI9A1) - Azide and BSA Free [NBP2-73935] - Staining of paraffin-embedded Human thyroid tissue using anti-RPA2 mouse monoclonal antibody.




Flow Cytometry: RPA2 Antibody (OTI9A1) - Azide and BSA Free [NBP2-73935] - Analysis of Jurkat cells, using anti-RPA2 antibody, (Red), compared to a nonspecific negative control antibody (Blue).



Western Blot: RPA2 Antibody (OTI9A1) - Azide and BSA Free [NBP2-73935] - Analysis of extracts (35ug) from 9 different cell lines by using g anti-RPA2 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).



Western Blot: RPA2 Antibody (OTI9A1) - Azide and BSA Free [NBP2-73935] - Ajuba associates with RPA in HTC75 and IMR90 cells. IP-Western probed with a total anti-RPA32 antibody after immunoprecipitation with a commercial anti-Ajuba antibody (com.) or anti-peptide serum from IMR90 or HTC75 extracts as indicated. Left panel: tota



Immunohistochemistry: RPA2 Antibody (OTI9A1) - Azide and BSA Free [NBP2-73935] - Staining of paraffin-embedded Adenocarcinoma of Human colon tissue using anti-RPA2 mouse monoclonal antibody.




Immunohistochemistry: RPA2 Antibody (OTI9A1) - Azide and BSA Free [NBP2-73935] - Staining of paraffin-embedded Human colon tissue using anti-RPA2 mouse monoclonal antibody.




Immunohistochemistry: RPA2 Antibody (OTI9A1) - Azide and BSA Free [NBP2-73935] - Staining of paraffin-embedded Human prostate tissue using anti-RPA2 mouse monoclonal antibody.



Flow Cytometry: RPA2 Antibody (OTI9A1) - Azide and BSA Free [NBP2-73935] - Analysis of Hela cells, using anti-RPA2 antibody, (Red), compared to a nonspecific negative control antibody (Blue).



Immunoprecipitation: RPA2 Antibody (OTI9A1) - Azide and BSA Free [NBP2-73935] - (Negative control: IP without adding anti-RPA2 antibody.). For each experiment, 500ul of DDK tagged RPA2 overexpression lysates (at 1:5 dilution with HEK293T lysate), 2ug of anti-RPA2 antibody and 20ul (0.1mg) of goat anti-mouse conjugated magnetic beads were mixed and incubated overnight. After extensive wash to remove any non-specific binding, the immuno-precipitated products were analyzed with rabbit anti-DDK polyclonal antibody.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-73935**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-73935](http://www.novusbio.com/reviews/submit/NBP2-73935)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

