

# Product Datasheet

## TSC22/TSC22D1 Antibody (OTI4H8) [PE] NBP2-74654PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-74654PE](http://www.novusbio.com/NBP2-74654PE)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-74654PE](http://www.novusbio.com/reviews/destination/NBP2-74654PE)



**NBP2-74654PE**

TSC22/TSC22D1 Antibody (OTI4H8) [PE]

**Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 ml                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C in the dark.                                                                    |
| <b>Clonality</b>     | Monoclonal                                                                                  |
| <b>Clone</b>         | OTI4H8                                                                                      |
| <b>Preservative</b>  | 0.05% Sodium Azide                                                                          |
| <b>Isotype</b>       | IgG2a                                                                                       |
| <b>Conjugate</b>     | PE                                                                                          |
| <b>Purity</b>        | Immunogen affinity purified                                                                 |
| <b>Buffer</b>        | PBS                                                                                         |

**Product Description**

|                    |                                                                                               |
|--------------------|-----------------------------------------------------------------------------------------------|
| <b>Host</b>        | Mouse                                                                                         |
| <b>Gene ID</b>     | 8848                                                                                          |
| <b>Gene Symbol</b> | TSC22D1                                                                                       |
| <b>Species</b>     | Human, Mouse, Rat                                                                             |
| <b>Immunogen</b>   | Full length recombinant protein of human TSC22/TSC22D1 (NP_006013) produced in HEK293T cells. |

**Product Application Details**

|                              |                                                                        |
|------------------------------|------------------------------------------------------------------------|
| <b>Applications</b>          | Western Blot, Flow Cytometry, Immunohistochemistry, CyTOF-ready        |
| <b>Recommended Dilutions</b> | Western Blot, Flow Cytometry, Immunohistochemistry, CyTOF-ready        |
| <b>Application Notes</b>     | Optimal dilution of this antibody should be experimentally determined. |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP2-74654PE**

---

|                    |                                                      |
|--------------------|------------------------------------------------------|
| H00008848-P01-10ug | Recombinant Human TSC22/TSC22D1 GST (N-Term) Protein |
| DKK300             | Kallikrein 3/PSA [HRP]                               |
| H00008848-P02-10ug | Recombinant Human TSC22/TSC22D1 GST (N-Term) Protein |
| AF3797             | Smad2/3 Antibody [Unconjugated]                      |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-74654PE](http://www.novusbio.com/reviews/submit/NBP2-74654PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

