

Product Datasheet

UBOX5 Antibody (OT11C5) - Azide and BSA Free NBP2-74760

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-74760

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-74760



NBP2-74760

UBOX5 Antibody (OT11C5) - Azide and BSA Free


Product Information	
Unit Size	100 ug
Concentration	LYOPH mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OT11C5
Preservative	No Preservative
Reconstitution Instructions	we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Purity	Immunogen affinity purified
Buffer	Lyophilized from PBS (pH 7.3) with 8% Trehalose

Product Description	
Description	Novus Biologicals Mouse UBOX5 Antibody (OT11C5) - Azide and BSA Free (NBP2-01077) is a monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	22888
Gene Symbol	UBOX5
Species	Human
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-130 and C 419-487 of human UBOX5(NP_955447) produced in E.coli.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, CyTOF-ready
Recommended Dilutions	Western Blot 1:200-500, Flow Cytometry 1:100, Immunohistochemistry 1:150, Immunocytochemistry/ Immunofluorescence 1:100, Immunohistochemistry-Paraffin, CyTOF-ready

Images


Western Blot: UBOX5 Antibody (OT11C5) - Azide and BSA Free [NBP2-74760] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY UBOX5 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UBOX5.




Immunocytochemistry/Immunofluorescence: UBOX5 Antibody (OT11C5) - Azide and BSA Free [NBP2-74760] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY UBOX5 .




Immunohistochemistry: UBOX5 Antibody (OT11C5) - Azide and BSA Free [NBP2-74760] - Staining of paraffin-embedded Carcinoma of Human prostate tissue using anti-UBOX5 mouse monoclonal antibody.




Flow Cytometry: UBOX5 Antibody (OT11C5) - Azide and BSA Free [NBP2-74760] - HEK293T cells transfected with either overexpression plasmid (Red) or empty vector control plasmid (Blue) were immunostaining by anti-UBOX5 antibody, and then analyzed by flow cytometry.



Immunohistochemistry: UBOX5 Antibody (OT11C5) - Azide and BSA Free [NBP2-74760] - Staining of paraffin-embedded Adenocarcinoma of Human endometrium tissue using anti-UBOX5 mouse monoclonal antibody.



Immunohistochemistry: UBOX5 Antibody (OT11C5) - Azide and BSA Free [NBP2-74760] - Staining of paraffin-embedded Carcinoma of Human bladder tissue using anti-UBOX5 mouse monoclonal antibody.




Immunohistochemistry: UBOX5 Antibody (OT11C5) - Azide and BSA Free [NBP2-74760] - Staining of paraffin-embedded Carcinoma of Human liver tissue using anti-UBOX5 mouse monoclonal antibody.




Immunohistochemistry: UBOX5 Antibody (OT11C5) - Azide and BSA Free [NBP2-74760] - Staining of paraffin-embedded Human endometrium tissue using anti-UBOX5 mouse monoclonal antibody.



Flow Cytometry: UBOX5 Antibody (OT11C5) - Azide and BSA Free [NBP2-74760] - Analysis of Hela cells, using anti-UBOX5 antibody, (Red), compared to a nonspecific negative control antibody (Blue).



Flow Cytometry: UBOX5 Antibody (OT11C5) - Azide and BSA Free [NBP2-74760] - Analysis of Jurkat cells, using anti-UBOX5 antibody, (Red), compared to a nonspecific negative control antibody (Blue).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-74760

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-74760

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

