

Product Datasheet

USP7 Antibody (OT11F12) [mFluor Violet 610 SE] NBP2-74814MFV610

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-74814MFV610

Updated 11/9/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-74814MFV610



NBP2-74814MFV610

USP7 Antibody (OT11F12) [mFluor Violet 610 SE]

| Product Information | |
|---------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | OT11F12 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 |
| Conjugate | mFluor Violet 610 SE |
| Purity | Immunogen affinity purified |
| Buffer | 50mM Sodium Borate |

| Product Description | |
|---------------------|---|
| Host | Mouse |
| Gene ID | 7874 |
| Gene Symbol | USP7 |
| Species | Human, Mouse, Rat |
| Reactivity Notes | Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions. |
| Immunogen | Full length human recombinant protein of human USP7(NP_003461) produced in HEK293T cell. |

| Product Application Details | |
|-----------------------------|--|
| Applications | Western Blot, Flow Cytometry, Immunohistochemistry |
| Recommended Dilutions | Western Blot, Flow Cytometry, Immunohistochemistry |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |

Images

USP7 Antibody (OT11F12) [mFluor Violet 610 SE] [NBP2-74814MFV610]
 - Vial of mFluor Violet 610 conjugated antibody. mFluor Violet 610 is optimally excited at 421 nm by the Violet laser (405 nm) and has an emission maximum of 613 nm.



mFluor™ Violet 610

| | |
|---------------------|-------------------|
| LASER (nm) | FILTER |
| Violet (405) | 605/30 |
| EXCITATION MAX (nm) | EMISSION MAX (nm) |
| 421 | 613 |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-74814MFV610

| | |
|---------------|---------------------------------------|
| NBP1-83293PEP | USP7 Recombinant Protein Antigen |
| E-519-025 | USP7 [Unconjugated] |
| NB100-56519 | DNMT1 Antibody (60B1220.1) - BSA Free |
| NB100-59787 | PML Protein Antibody - BSA Free |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-74814MFV610

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

