

Product Datasheet

Galectin-8 Antibody (JB85-35)

NBP2-75501

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-75501

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-75501



NBP2-75501**Galectin-8 Antibody (JB85-35)**


Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	JB85-35
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	TBS (pH7.4), 0.05% BSA, 40% Glycerol
Target Molecular Weight	40 kDa

Product Description	
Description	Novus Biologicals Rabbit Galectin-8 Antibody (JB85-35) (NBP2-75501) is a recombinant monoclonal antibody validated for use in IHC, WB, Flow, ICC/IF and IP. Anti-Galectin-8 Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3964
Gene Symbol	LGALS8
Species	Human, Mouse, Rat
Immunogen	Synthetic peptide within Human Galectin-8 aa 1-50 / 317. (SwissProt: O00214 Human; SwissProt: Q9JL15 Mouse; SwissProt: Q62665 Rat)


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:1000, Flow Cytometry 1:50-1:100, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1:400-1:800, Immunoprecipitation, Immunohistochemistry-Paraffin 1:50-1:200

Images


Western Blot: Galectin-8 Antibody (JB85-35) [NBP2-75501] - Analysis of Galectin 8 on mouse testis and rat testis tissue lysates using anti-Galectin 8 antibody at 1/500 dilution.




Immunocytochemistry/Immunofluorescence: Galectin-8 Antibody (JB85-35) [NBP2-75501] - Staining Galectin 8 in HepG2 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Immunohistochemistry-Paraffin: Galectin-8 Antibody (JB85-35) [NBP2-75501] - Analysis of paraffin-embedded human colon cancer tissue using anti-Galectin 8 antibody. Counter stained with hematoxylin.




Flow Cytometry: Galectin-8 Antibody (JB85-35) [NBP2-75501] - Analysis of PC-3M cells with Galectin 8 antibody at 1/100 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti-rabbit IgG was used as the secondary antibody.




Immunocytochemistry/Immunofluorescence: Galectin-8 Antibody (JB85-35) [NBP2-75501] - Staining Galectin 8 in A431 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



Immunohistochemistry-Paraffin: Galectin-8 Antibody (JB85-35) [NBP2-75501] - Analysis of paraffin-embedded human kidney tissue using anti-Galectin 8 antibody. Counter stained with hematoxylin.



Immunohistochemistry-Paraffin: Galectin-8 Antibody (JB85-35) [NBP2-75501] - Analysis of paraffin-embedded human liver tissue using anti-Galectin 8 antibody. Counter stained with hematoxylin.



Publications

Jia X, Knyazeva A, Zhang Y et al. V. cholerae MakA is a cholesterol-binding pore-forming toxin that induces non-canonical autophagy The Journal of cell biology 2022-12-05 [PMID: 36194176] (ICC/IF, Human)

Details:

Dilution used in ICC/IF 1:400

Roberts JA, Varma VR, An Y Et al. A brain proteomic signature of incipient Alzheimer's disease in young APOE epsilon4 carriers identifies novel drug targets Science advances 2021-11-12 [PMID: 34757788] (WB, Human, Mouse)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-75501

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-55285PEP	Galectin-8 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-75501

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

